NRC Publications Archive Archives des publications du CNRC

COVID-Net USPro: an explainable few-shot deep prototypical network for COVID-19 screening using point-of-care ultrasound

Song, Jessy; Ebadi, Ashkan; Florea, Adrian; Xi, Pengcheng; Tremblay, Stéphane; Wong, Alexander

This publication could be one of several versions: author's original, accepted manuscript or the publisher's version. / La version de cette publication peut être l'une des suivantes : la version prépublication de l'auteur, la version acceptée du manuscrit ou la version de l'éditeur.

For the publisher's version, please access the DOI link below./ Pour consulter la version de l'éditeur, utilisez le lien DOI ci-dessous.

Publisher's version / Version de l'éditeur:

https://doi.org/10.3390/s23052621

Sensors, 23, 5, pp. 1-18, 2023-02-27

NRC Publications Archive Record / Notice des Archives des publications du CNRC : https://publications-cnrc.canada.ca/fra/voir/objet/?id=abe11456-3345-4f98-85d0-486109e8a0cd

Access and use of this website and the material on it are subject to the Terms and Conditions set forth at https://nrc-publications.canada.ca/eng/copyright

READ THESE TERMS AND CONDITIONS CAREFULLY BEFORE USING THIS WEBSITE.

L'accès à ce site Web et l'utilisation de son contenu sont assujettis aux conditions présentées dans le site https://publications-cnrc.canada.ca/fra/droits

LISEZ CES CONDITIONS ATTENTIVEMENT AVANT D'UTILISER CE SITE WEB.

Questions? Contact the NRC Publications Archive team at

PublicationsArchive-ArchivesPublications@nrc-cnrc.gc.ca. If you wish to email the authors directly, please see the first page of the publication for their contact information.

Vous avez des questions? Nous pouvons vous aider. Pour communiquer directement avec un auteur, consultez la première page de la revue dans laquelle son article a été publié afin de trouver ses coordonnées. Si vous n'arrivez pas à les repérer, communiquez avec nous à PublicationsArchive-ArchivesPublications@nrc-cnrc.gc.ca.







MDPI

Article

COVID-Net USPro: An Explainable Few-Shot Deep Prototypical Network for COVID-19 Screening Using Point-of-Care Ultrasound

Jessy Song ¹, Ashkan Ebadi ^{1,2,*}, Adrian Florea ³, Pengcheng Xi ^{1,4}, Stéphane Tremblay ⁴ and Alexander Wong ^{1,5}

- Department of Systems Design Engineering, University of Waterloo, Waterloo, ON N2L 3G1, Canada
- ² Digital Technologies Research Centre, National Research Council Canada, Toronto, ON M5T 3J1, Canada
- Department of Emergency Medicine, McGill University, Montreal, QC H4A 3J1, Canada
- ⁴ Digital Technologies Research Centre, National Research Council Canada, Ottawa, ON K1A 0R6, Canada
- ⁵ Waterloo Artificial Intelligence Institute, Waterloo, ON N2L 3G1, Canada
- * Correspondence: ashkan.ebadi@nrc-cnrc.gc.ca

Abstract: As the Coronavirus Disease 2019 (COVID-19) continues to impact many aspects of life and the global healthcare systems, the adoption of rapid and effective screening methods to prevent the further spread of the virus and lessen the burden on healthcare providers is a necessity. As a cheap and widely accessible medical image modality, point-of-care ultrasound (POCUS) imaging allows radiologists to identify symptoms and assess severity through visual inspection of the chest ultrasound images. Combined with the recent advancements in computer science, applications of deep learning techniques in medical image analysis have shown promising results, demonstrating that artificial intelligence-based solutions can accelerate the diagnosis of COVID-19 and lower the burden on healthcare professionals. However, the lack of large, well annotated datasets poses a challenge in developing effective deep neural networks, especially in the case of rare diseases and new pandemics. To address this issue, we present COVID-Net USPro, an explainable few-shot deep prototypical network that is designed to detect COVID-19 cases from very few ultrasound images. Through intensive quantitative and qualitative assessments, the network not only demonstrates high performance in identifying COVID-19 positive cases, using an explainability component, but it is also shown that the network makes decisions based on the actual representative patterns of the disease. Specifically, COVID-Net USPro achieves 99.55% overall accuracy, 99.93% recall, and 99.83% precision for COVID-19-positive cases when trained with only five shots. In addition to the quantitative performance assessment, our contributing clinician with extensive experience in POCUS interpretation verified the analytic pipeline and results, ensuring that the network's decisions are based on clinically relevant image patterns integral to COVID-19 diagnosis. We believe that network explainability and clinical validation are integral components for the successful adoption of deep learning in the medical field. As part of the COVID-Net initiative, and to promote reproducibility and foster further innovation, the network is open-sourced and available to the public.

Keywords: ultrasonic imaging; lung; COVID-19; few-shot learning; deep explainable architecture



Citation: Song, J.; Ebadi, A.; Florea, A.; Xi, P.; Tremblay, S.; Wong, A. COVID-Net USPro: An Explainable Few-Shot Deep Prototypical Network for COVID-19 Screening Using Point-of-Care Ultrasound. *Sensors* 2023, 23, 2621. https://doi.org/10.3390/s23052621

Academic Editor: Yu-Dong Zhang

Received: 26 January 2023 Revised: 20 February 2023 Accepted: 21 February 2023 Published: 27 February 2023



Copyright: © 2023 by the authors. Licensee MDPI, Basel, Switzerland. This article is an open access article distributed under the terms and conditions of the Creative Commons Attribution (CC BY) license (https://creativecommons.org/licenses/by/4.0/).

1. Introduction

The Coronavirus Disease 2019, or COVID-19, caused by severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2), has been continuously impacting individuals' wellbeing and the global healthcare systems [1]. Despite the vaccination efforts, policies, and regulations in place, due to the rapid transmission of the virus and waves of rising cases, the development of effective screening and risk stratification methods remains to be a critical need in controlling the disease [2]. Various types of diagnostic tools, including reverse transcription-polymerase chain reaction (RT-PCR), rapid antigen detection tests,

Sensors **2023**, 23, 2621 2 of 18

and antibody tests, have been developed and adapted globally to increase the rate of screening. While RT-PCR has been the gold standard test for diagnosing COVID-19, the technique involves large labour and laboratory resources and is time-consuming [3]. Other rapid antigen tests and antibody tests with varying sensitivity are also less reliable in comparison to RT-PCR tests [3].

For people with significant respiratory symptoms, medical imaging is used to identify the disease and assess the severity of the disease progression [4]. Under this protocol, a computed tomography (CT) scan, chest X-ray (CXR), or point-of-care ultrasound (POCUS) imaging can be performed and used clinically as an alternative diagnostic tool [2]. To make a diagnosis, acute care physicians and radiologists visually inspect the radiographic images to find patterns related to symptoms and to assess the severity of COVID-19 infection and deformation [3]. During times of high transmission rate of COVID-19, the large influx of patients increases the burden on clinicians and radiologists. Medical image processing and artificial intelligence (AI) can assist in reducing this burden and accelerate the diagnostic and decision-making process, as existing models and algorithms continue to improve and the amount of available medical image data continues to grow [5–7].

Different imaging modalities, including CT scan, X-ray, and ultrasound, may be used in the diagnosis of COVID-19 and offer varying diagnostic values [8]. Chest CT scan is the most sensitive imaging modality in the initial diagnosis and management of confirmed cases, but it is more expensive and time-consuming [5,8]. In contrast, ultrasound imaging is more accessible and portable, cheap, and safer, as radiation is not involved during the examination, which are desirable properties for its usage, especially in resource-limited settings/environments/areas/regions [8].

Deep learning usually requires a large set of training examples [4,7,9]. However, due to the nature of novel diseases, the availability of such a huge amount of well annotated data poses a great challenge to learning algorithms. Few-shot learning is an approach where a model is trained to classify new data based on a limited number of samples exposed in training [10]. This approach resembles how humans learn, as we can recognize new object classes from very few instances, and it is different from conventional deep neural networks that require a large amount of data in the training phase [10]. Since the few-shot model requires much fewer data to train, the computational costs are also significantly reduced [10]. These properties make it an appropriate and promising approach for COVID-19 as it relates to novel and rare disease diagnosis. One approach for few-shot learning is metric-based learning. As a few-shot metric-based learning approach, prototypical networks (PN) perform classification by computing distances to prototype representations of each class [10]. PN has shown state-of-the-art (SOTA) results on other datasets and domains (e.g., [11–13]), proving that simple design decisions can yield significant improvements over other complicated architectures and meta-learning approaches [10].

In this work, we present an open-source explainable deep prototypical network, called COVID-Net USPro, that can detect COVID-19 cases with high accuracy, precision, and recall from a very limited number of lung ultrasound (LUS) images. When trained with only 5 shots, COVID-Net USPro classifies between positive and negative COVID-19 cases with 99.55% overall accuracy, 99.93% recall, and 99.83% precision for COVID-19 positive cases. Intensive experimentation was conducted (e.g., testing different image encoders, varying training conditions, and the number of classes to optimize the network) to assess the performance of COVID-Net USPro. To ensure the network's fairness and accountability, an explainability module is constructed to assess the network's decisions with visual explanation tools, i.e., Grad-CAM [14] and GSInquire [15]. Moreover, our contributing clinician (A.F.) carefully verified and validated the pipeline, as well as the results, to ensure the validity of the proposed solution from the clinical perspective. To facilitate the adoption and openness of AI in healthcare, support reproducibility, and encourage innovation, the network, and all the experiment scripts, are open-sourced at the project's Github repository.

Our work contributes to the existing body of literature on rare disease medical image analysis using few-shot learning in at least the three following ways:

Sensors **2023**, 23, 2621 3 of 18

1. We present a few-shot network that reaches 99.55% accuracy when trained with only five shots, while other related works achieving similar (or lower) performance require larger numbers of training examples (see details in Section 2).

- 2. COVID-Net USPro is an explainable network, as demonstrated quantitatively by analysis from two explainability visualization tools and qualitatively by our clinician validation.
- COVID-Net USPro is open-sourced and available to the public, which helps promote
 the reproducibility and accessibility of AI in healthcare. This would encourage further
 innovation in the field of deep learning applied to medical image analysis for novel
 disease diagnosis.

The remainder of this paper is structured as follows. Section 2 describes related previous studies and highlights their limitations. Section 3 explains data, techniques, and the experiments conducted to assess the network performance in detail. Section 4 presents findings of the research, including quantitative performance assessment and qualitative explainability analysis results. Results are then discussed and summarized in Section 5, where some limitations of the research and future directions are also presented.

2. Related Work

There are several studies that aim to apply deep learning to the screening and detection of COVID-19 cases. As an open-source and open-access initiative, the COVID-Net [5,7,9,16] includes research on the application of deep learning neural networks using a multitude of image modalities, such as CT, X-ray, and ultrasound images. Multiple works have demonstrated the effectiveness of deep learning in the classification of CT and X-ray images. For example, Aboutalebi et al. proposed the COVID-Net CXR network [17], which is a tailored deep convolutional neural network (DCNN/CNN) for the detection of COVID-19 cases. The network that was trained on 5210 chest X-ray images achieved an overall accuracy of 98.3% and 97.5% sensitivity for COVID-19 cases. In another work, Ozturk et al. [18] proposed a DCNN network that was previously used for the you-only-look-once (YOLO) real-time object detection system to classify X-ray images. Their proposed network achieved 98.08% accuracy for binary COVID-19 case detection. Afshar et al. [19] proposed a capsule network, called COVID-CAPS, which achieved over 98% accuracy and, specificity, trained on 112,120 X-ray images. Gunraj et al. [6] proposed the COVID-Net CT network for COVID-19 detection from CT images. Being trained on 194,922 CT images, the network scored 96.2% in sensitivity and 99% in specificity for COVID-19 cases. The potential of including both CT-scan and X-ray images for classification is also explored. For instance, Thakur [20] presents a DCNN-based model achieving over 99% accuracy and precision for COVID-19 detection using 11,095 X-ray and CT images. Few works proposed DL-based solutions for COVID-19 detection from ultrasound images. In a recent study, MacLean et al. [7] proposed a highly efficient self-attention deep neural network, called COVID-Net US, for COVID-19 detection from POCUS ultrasound imaging. The network was trained on 3947 positive COVID-19 and 3697 negative normal case images and achieved an area under the receiver operating curve (AUC) of over 98%. In another study, Diaz-Escobar et al. [21] leveraged pre-trained neural networks, such as VGG19 [22], InceptionV3 [23], and ResNet50 [24], in the detection of COVID-19 from ultrasound images and achieved 89.1% accuracy and AUC of 97.1%. One main requirement of conventional deep neural network architectures in most of the existing research is a large amount of training data, as in all the mentioned works above, datasets surpassed 5000 total images [7,9,17]. This may pose a serious limitation in the case of novel/rare diseases and new pandemics where not many images are available for training the DL model.

The application of few-shot learning techniques has also been investigated. For example, Shorfuzzaman et al. [25] proposed MetaCOVID, a Siamese few-shot neural network with the contrastive loss for detecting COVID-19 using CXR images. The performance of the best network achieved an accuracy of 95.6% and an AUC of 97% when tested in a 10-shot setting in the inference phase. In another work by Ebadi et al. [26], a deep Siamese convolutional network called COVID-Net FewSE was proposed that can detect COVID-19

Sensors **2023**, 23, 2621 4 of 18

positive cases with 90% recall and accuracy of 99.93% when the network is provided with only 50 examples in the training phase. In the work by Karnes et al. [27], the use of adaptive few-shot learning for ultrasound COVID-19 detection was studied.

Although the feasibility of adopting few-shot learning techniques for COVID-19 detection from medical imaging has been already investigated, most of the previous studies only focused on the performance of the network, analyzing its accuracy, precision, and recall. A comprehensive quantitative analysis of network explainability to ensure that decisions are made based on actual patterns is either missing or inadequate. In addition, a thorough validation assessment, performed by clinicians as the domain experts, is also lacking. There are limited discussions on whether the data interpretation process of these well performing networks aligns with real clinical settings [25]. These limitations could jeopardize the full understanding of the network and hinders the adoption of the network in the real clinical environment. In this work, we aim to address these limitations by presenting a high-performing network that can accurately detect COVID-19 cases when presented with only five shots. We also perform a comprehensive explainability analysis to validate network behaviour. Additionally, most importantly, we further validate the network and findings by an experienced practicing clinician.

3. Materials and Methods

3.1. Data

We use the COVIDx-US dataset v1.4. [1] as the data source. COVIDx-US is an openaccess benchmark dataset of lung ultrasound imaging data that contains 242 videos and 29,651 processed images of patients with COVID-19 infection, non-COVID-19 infection (mainly pneumonia), other lung conditions, and normal control cases. The dataset provides LUS images captured with two kinds of probes, linear probes, which produce a square or rectangular image, and convex probes, which allow for a wider field of view [28]. Due to the difference in the field of view, combining the linear and convex probe data in training may increase noise and influence the performance of the network. As there are also a low number of COVID-19 positive examples captured with the linear probes in the dataset, we exclude them from this study. A total number of 25,262 convex LUS images are then randomly split into the training set, containing 90% of images, and the unseen test set with the remaining 10% of images, ensuring all frames from each video are either in train or test set to avoid data leakage. The training set is then split into 80–20%, representing the training and validation datasets. The validation dataset is used for hyperparameter tuning and performance assessment in the training phase. All images are re-scaled to 224×224 pixels to keep the images across the entire dataset consistent. The dataset is further augmented by rotating each image by 90°, 180°, and 270°, resulting in a total of 101,048 images (25, 262×4). This rotation technique is an appropriate method for increasing the dataset size, as it keeps the images and areas of interest for clinical decisions unaltered and in-bound [29].

3.2. *Methods*

COVID-Net USPro is a prototypical few-shot learning network that trains in an episodic learning setting. It uses a distance metric for assessing similarities between a set of unlabelled data, i.e., query set, and labelled data, i.e., support set. Labelled data can be used to compute a single *prototype* representation of the class, and unlabelled data are assigned to the class of the prototype they are closest to. A prototypical network [10] is based on the idea that there exists an embedding in which points in a class cluster around a single prototype representation for the class. During the training phase, a neural network is used to learn the non-linear mapping of the inputs to an embedding space, and a class prototype is computed as the mean of its support set data in the embedding space. Classification is then done by finding the nearest class prototype for each query point based on a specified distance metric. An episodic approach is used to train the model, where in each training episode, the few-shot task is simulated by sampling the data point in mini-batches

Sensors **2023**, 23, 2621 5 of 18

to make the training process consistent with the testing environment. The performance of the network is evaluated using the unseen test set. Both quantitative performance analysis based on accuracy, precision, and recall and qualitative explainability analysis are conducted. A high-level conceptual flow of the analysis is presented in Figure 1.




Figure 1. High-level conceptual flow of the analysis.

We define the classification problem as a K-way N-shot episodic task, where K denotes the number of classes present in the dataset, and N denotes the number of images for each class in each episode. For a given dataset, N images from each of the K classes are sampled to form the support set, and another M images from each class are sampled to form the query set. The network then aims to classify the images of the query set based on the K*N total images presented in the support set. In this work, we consider three classification scenarios and formulate the problem as 2-way, 3-way, and 4-way classification problems. Details are included under Section 3.3.3.

The few-shot classification with a prototypical network can be summarized into three steps: (1) encoding of the images, (2) generating class prototypes, and (3) assigning labels to query samples based on distances to the class prototypes. Let $S = \{(\mathbf{x}_{(1,s)}, y_{(1,s)}), \ldots, (\mathbf{x}_{(N,s)}, y_{(N,s)})\}$ and $Q = \{(\mathbf{x}_{(1,q)}, y_{(1,q)}), \ldots, (\mathbf{x}_{(N,q)}, y_{(N,q)})\}$ be the support and query sets, respectively, where each $x_i \in R^D$ is a D-dimensional example feature vector and $y_i \in \{1, \ldots K\}$ is the label of the example. The prototypical network contains an image encoder $f_\phi: R^D \to R^H$ that transforms each image x_i onto a H-dimensional embedding space where images of the same class cluster together. Class prototypes are then generated for each class by averaging the embedding image vectors in the support set, where $v_k = \frac{1}{N} \sum_{i=1}^N f_\phi(x_{i,s}^{(k)})$ denotes the prototype of class k [10]. To classify the query image, a distance metric is used where distances between the embedding vector of a query image and each of the class prototypes are computed. In this work, the squared Euclidean distance $d(v,q) = ||v-q|| = \sqrt{\sum (v_i - q)^2}$

Sensors **2023**, 23, 2621 6 of 18

is used, where q is the embedding vector of the query image, and v_i is the embedding vector of the i-th prototype. The choice of the squared Euclidean distance instead of other distance metrics, e.g., cosine distance, is validated by Snell et al. [10], who proved that metrics that are Bregman divergences, e.g., euclidean distance metrics, perform better in the calculation of class prototypes based on embeddings in prototypical networks. After distances are computed, a SoftMax function is applied over the distances to the prototypes to compute the probabilities of the query image being in each class. The class with the highest probability is then assigned to the query image.

In the training phase, the network learns by minimizing a loss function, i.e., the negative log-SoftMax function ($J = -\log{(p(y=k|x))}$) of the true class k. An Adam optimizer with an initial learning rate of 0.001 is used, and the learning rate is reduced if the validation loss is not improved after 3 epochs. In each episode, a subset of data points is randomly selected, forming a support and query set. The loss terms on training and validation sets are calculated at the end of each training episode. To facilitate an effective training process and prevent over-fitting, early stopping is implemented to stop the training process when the validation loss is not improved after 5 epochs. A total of 10 epochs is set for all training processes, and 200 episodes are set for each training epoch. Figure 2 presents an architecture design overview of the COVID-Net USPro network.




Figure 2. COVID-Net USPro, network architecture design.

The trained model's performance is evaluated quantitatively and qualitatively. In the quantitative analysis, the model's accuracy, precision, and recall for each class are analyzed and reported. In the qualitative analysis, model explainability is investigated and visualized. Explainable artificial intelligence (XAI) has been an important criterion when assessing whether neural networks can be applied to clinical settings [30]. While AI-driven systems may show high accuracy and precision in analyzing medical images, lack of reasonable explainability will spark criticism of the network's adoption [30]. COVID-Net USPro's explainability is assessed using two approaches, i.e., gradient-weighted class activation map (Grad-CAM) [14] and GSInquire [15], on a selected dataset containing correctly classified COVID-19 and normal cases with high confidence (i.e., >99.9% probability), as well as falsely predicted COVID-19 and normal cases. Grad-CAM generates a visual explanation of the input image using the gradient information flowing into the last convolutional layer of the convolutional neural network (CNN) encoder and assigns importance values to each neuron for making a classification decision [14]. The output is a heatmap-overlayed image that shows the image regions that impact the particular classification decision made by the network [14]. The other tool, GSInquire, identifies the critical factors in an input image that are shown to be integral to the decisions made by the network in a generative synthesis approach [15]. The result is an annotated image highlighting the critical region, which could drastically change the classification result if removed [15]. Results from both tools are reviewed by our contributing clinician (A.F.) with experience in ultrasound image analysis to assess whether clinically important patterns are captured by the network.

Sensors **2023**, 23, 2621 7 of 18

3.3. Experiment Settings

We comprehensively assess the performance of COVID-Net USPro in detecting COVID-19 cases from ultrasound images by testing various training conditions, such as different image encoders, the number of shots available for training, and classification task types. Details are further discussed in this section.

3.3.1. Image Encoders

To leverage the power of transfer learning, we experiment with multiple encoders, including, but not limited to, the ResNet, DenseNet, and VGG networks [22,24,31]. Pretrained models refer to using model parameters pre-trained on ImageNet [32]. To concisely summarize the main results, we report the top-4 performing encoders with respect to our research objectives:

- ResNet18L1: Pre-trained ResNet18 [24], with trainable parameters on the final connected layer and setting out features as the number of classes. We consider this pre-trained network as the baseline model for encoders, as it contains the least number of layers and retrained parameters.
- ResNet18L5: Pre-trained ResNet18 [24], with trainable parameters on the last four convolutional layers and final connected layer. Out features set to the number of classes.
- **ResNet50L1:** Pre-trained ResNet50 [24], with trainable parameters on the final connected layer and setting out features as the number of classes.
- ResNet50L4: Pre-trained ResNet50 [24], with trainable parameters on the last three
 convolutional layers and final connected layer. Out features set to the number
 of classes.

3.3.2. Number of Training Shots

The optimal number of shots for maximized performance is tested by training the model under 5, 10, 20, 30, 40, 50, 75, and 100-shot scenarios. For selected models showing a steady increase of performance over increasing shots, 150 and 200-shot conditions are also tested to further verify that the maximum performance is reached at 100-shot. To ensure the training process is faithful to the testing environment, the number of shots for each class presented in each episode is the same in the support and query sets in both training and test phases. For example, in the 5-shot scenario, five images in each class are presented for both the support set and the query set in the training phase, and the same follows in the test phase.

3.3.3. Problem Formulation

In comparison to other classes, e.g., non-COVID-19 and normal cases, the ability of the model to correctly identify COVID-19-positive cases is valued the most. The classification problem for identifying COVID-19 is formulated in three different scenarios as follows, in ascending order of data complexity:

- **2-way classification:** Data from all three other classes, namely, the 'normal' class, 'non-COVID-19' class, and 'other' class, are viewed as a combined COVID-19 negative class. The network learns from COVID-19 positive and COVID-19 negative datasets in this setting.
- **3-way classification:** As the 'other' class contains data from multiple different lung conditions, it has the highest variations and may disrupt the network's learning process due to the lack of uniformity in the data. In the 3-class classification, the 'other' class is excluded, and the network is trained to classify the remaining three classes.
- **4-way classification:** As the dataset contains four classes, the four-class classification condition remains in this setting, and the network is trained to classify 'COVID-19', 'normal', 'non-COVID-19', and 'other' classes.

The network hyperparameters and training settings are listed in Appendix A.

Sensors **2023**, 23, 2621 8 of 18

4. Results

This section summarizes the quantitative performance results of all combinations of experiment settings listed in Section 3.3, as well as the results of the network explainability analysis.

4.1. Quantitative Performance Analysis

The performance of COVID-Net USPro is evaluated using the overall accuracy, precision, and recall for each class. As the performance of the model to detect COVID-19 cases is the most important for current clinical use cases, the following report contains precision and recall for the COVID-19 class only. To reduce table size, Table 1 only summarizes the performance of the network under 5-shot and 100-shot scenarios for encoders that scored over 80% across all evaluated metrics. The complete experiment results for all classes, encoders, and shot numbers are detailed in Appendix B.

Table 1. Summary of classification results for 5-shot and 100-shot conditions.

Scenario No. Shots		Model	Accuracy	Precision	Recall
2-way	5	ResNet18L1	0.9310	0.9384	0.9380
2-way	5	ResNet18L5	0.9955	0.9983	0.9993
2-way	5	ResNet50L1	0.9605	0.9567	0.9740
2-way	5	ResNet50L4	0.9910	0.9933	0.9900
2-way	100	ResNet18L1	0.9741	0.9741	0.9744
2-way	100	ResNet18L5	0.9999	1.0000	0.9999
2-way	100	ResNet50L1	0.9953	0.9954	0.9953
2-way	100	ResNet50L4	0.9999	0.9999	1.0000
3-way	5	ResNet18L1	0.9487	0.9470	0.9370
3-way	5	ResNet18L5	0.9947	0.9992	0.9920
3-way	5	ResNet50L1	0.9690	0.9761	0.9610
3-way	5	ResNet50L4	0.9953	0.9940	0.9940
3-way	100	ResNet18L1	0.9863	0.9858	0.9857
3-way	100	ResNet18L5	0.9998	0.9994	0.9998
3-way	100	ResNet50L1	0.9973	0.9947	0.9980
3-way	100	ResNet50L4	0.9998	0.9995	1.0000
4-way	5	ResNet18L1	0.8733	0.9283	0.898
4-way	5	ResNet18L5	0.9855	0.9980	0.9970
4-way	5	ResNet50L1	0.9320	0.9622	0.9550
4-way	5	ResNet50L4	0.9777	0.9967	0.9910
4-way	100	ResNet18L1	0.9369	0.9735	0.9643
4-way	100	ResNet18L5	0.9884	1.0000	1.0000
4-way	100	ResNet50L1	0.9811	0.9952	0.9962
4-way	100	ResNet50L4	0.9899	1.0000	1.0000

Across all classification scenarios and models, performance is higher under the 100-shot training condition than in the 5-shot condition, with performance metrics increasing from 5-shot and plateauing after 75-shot, as shown in Figure 3. The ResNet-based networks demonstrate the ability to classify COVID-19 with precision and recall above 87% consistently under both 5-shot, and above 99% under 100-shot conditions. As seen in Table 1, increasing the number of classes (in 3-way and 4-way classification scenarios) reduces the performance of the network. This is expected as the classification problem becomes more complex with more numbers of classes. However, this performance difference among the three classification scenarios is reduced when the number of shots increases, as more examples available for training improve the network's ability to distinguish between multiple classes. As shown in Figure 3, the performance of ResNet18L1, which has fewer layers than ResNet50L4, is lower than ResNet50L4 at early training. ResNet50L4 is expected to adapt better to input images in early training epochs and extract deeper and more complex representations. In addition, models with more fine-tuned final convolutional layers

Sensors **2023**, 23, 2621 9 of 18

(i.e., ResNet18L5 and ResNet50L4) achieve higher accuracy, precision, and recall. Therefore, it can be said that, while using pre-trained parameters and simpler models, reducing the computational complexity and space, tailoring parameters on the final convolutional layers to the ultrasound images, and deepening image encoding can boost performance.




Figure 3. Performance results with increasing shots trained under 4-way condition: (a) Pre-trained ResNet18 with trainable parameters on the final connected layer and setting out features as the number of classes (ResNet18L1). (b) Pre-trained ResNet50 with trainable parameters on the last 3 convolutional layers and final connected layer (ResNet50L4).

In the 2-way and 3-way classification scenarios, it is observed that the precision and recall for classes other than COVID-19 are similar to the COVID-19 class. In the 4-way classification scenario, the precision and recall for the 'other' class are \approx 2–3% lower than those for 'non-COVID-19', 'normal', and 'COVID-19' classes. This is expected since the 'other' class covers various lung conditions/diseases that encompass a larger range of image features and variations. Overall, with precision and recall achieving similar magnitude in the 2-way, 3-way, and 4-way classifications, the network also demonstrates the ability to distinguish between multiple classes/diseases. In comparison to studies outlined in Section 2, the performance of COVID-Net USPro networks, which are tailored to ultrasound images with re-trained parameters, is improved. Accuracy of ResNet50L1 and ResNet50L4 exceeds 98% under a 4-way 5-shot setting while past works such as the MetaCOVID proposed by Shorfuzzaman et al. [25], which also applied a few-shot approach, achieved 95.6% accuracy under a 3-way, 10-shot setting. Additionally, the sensitivity of COVID-Net USPro for COVID-19 cases is also higher than previously reported deep networks that were trained on many images of other modalities such as X-ray or CT, where they scored 97.5% in the best-performing case [6].

4.2. Clinical Validation and Network Explainability Analysis

In addition to the intensive quantitative performance analysis, we clinically validated the network outputs to ensure that the network captures important patterns in the ultrasound images. For this purpose, our contributing clinician (A.F.) reviewed a randomly selected set of images and reported his findings and observations. Our contributing clinician (A.F.) is an Assistant Professor in the Department of Emergency Medicine and the ultrasound co-director for undergraduate medical students at McGill University. He is practicing Emergency Medicine full-time at Saint Mary's Hospital in Montreal.

Figure 4 presents two examples of COVID-19 positive ultrasound images, annotated by Grad-CAM and GSInquire. As seen, the annotated images contain the lung pleura region at the top of the image, while the second example (Figure 4b) also marks the bottom region with high importance. B-lines, or the light comet-tail artifacts extending from pleura

Sensors **2023**, 23, 2621 10 of 18

to the bottom of the image, and the presence of dark regions interspacing the B-lines at the bottom part of the image, correspond to signs of lung consolidation and are indicators of abnormality [33]. Hence, the visual annotations for the second example (Figure 4b) are more representative of disease-related patterns within the ultrasound image. Figure 4a is one of the examples where the model considers the rib as a structure of interest, which is not an abnormality. Hence, although the model correctly classified the specific image, the decision was made based on invalid clinical factors.




Figure 4. COVID-19 positive case examples correctly classified by COVID-Net USPro with high confidence: (a) an example of wrong decision factors. (b) an example of a decision made based on disease-related patterns.

We experimented with two strategies to solve the above-mentioned issue. First, as the current dataset includes images of different qualities, we excluded images with low quality to evaluate their impact on explainability. These low-quality images refer to those with insufficient image depth or the lack of representative features. A severity grade was introduced by COVIDx-US dataset v1.4, called lung ultrasound score (LUSS), which rates each ultrasound video on a scale of 0 to 3 [1]. A score of 0 corresponds to the presence of only normal features, and 3 corresponds to the presence of severe disease artifacts [1]. In the first attempt to improve the network further, images from videos with a score of 0, representing the normal class, and images from videos with scores of 2 and 3, representing the COVID-19 class, are used to train a binary classification model. By observing the annotated images, the network shows to focus more on the bottom regions of the images, though cases, where the network focuses on the top pleura region, are still present. The second strategy to further improve model explainability is to exclude regions above the pleura (i.e., soft tissue) of the images so that the network focuses on the disease-defining features, which are mostly present at the bottom of the images below the lung pleura. Our experiments confirm the effectiveness of this strategy, as the network shows to focus mainly on the bottom regions of the images. Hence, combining the first and second strategies, a binary model with LUSS score filtered and cropped images is trained. Figure 5 shows examples from the analysis after cropping and filtering images. As suggested from the annotated examples and confirmed by our contributing clinician (A.F.), clinically relevant artifacts such as B-lines and lung consolidation are clearly identified in COVID-19-positive images by COVID-Net USPro after implementing the two strategies.

Sensors 2023, 23, 2621 11 of 18




Figure 5. Four cropped COVID-19 positive examples predicted correctly with high confidence by COVID-Net USPro (**a–d**), while recognizing disease artifacts, e.g., extended B-lines.

5. Discussion

Deep neural network architectures have shown promising results in a wide range of tasks, including predictive and diagnostic tasks. However, such networks require a massive amount of labelled data to train, which is against the nature of new pandemics and novel diseases where there are no or very few data samples available, especially in the initial stages. Using a diverse complex benchmark dataset, i.e., COVIDx-US, we introduced the COVID-Net USPro network, tailored to detect COVID-19 infection with high accuracy from very few ultrasound images. The proposed deep prototypical network leverages deep pre-trained models with fine-tuned parameters on final layers to reduce computational complexity and achieve high classification performance when only 5 examples from each class are presented during the training phase. Accuracy, precision, and recall for the best

Sensors **2023**, 23, 2621 12 of 18

performing network are over 99%, which are comparable to or outperforming other existing works, even those that used large-scale datasets to train the models [7,27]. As mentioned in Section 2, MetaCOVID network, proposed by Shorfuzzaman et al. [25], which adapted the few-shot learning approach achieved 95.6% accuracy when tested under a 3-way 10-shot condition. COVID-Net FewSE, proposed by Ebadi et al. [26], achieved over 99% accuracy when trained under a 3-way 50-shot setting. In comparison to these networks, COVID-Net USPro is able to achieve over 98.5% accuracy in 2-way, 3-way, and 4-way classification problem settings by using only 5 shots for training. These properties are not only highly crucial for the control of the COVID-19 pandemic, but also for screening patients for new/rare diseases or pandemics as the proposed network can be easily tuned.

In contrast to the lack of interpretability analysis in the previous works, we also intensively assessed the explainability of the network and clinically validated its performance. Experimental results demonstrate that COVID-Net USPro can not only achieve high performance in terms of accuracy, precision, and recall, but also shows predictive behaviour that is consistent with clinical interpretation, as validated by our contributing clinician (A.F.). In addition, as part of the explainability-driven performance validation process, we proposed and implemented two strategies to further improve the network performance in accordance with the background clinical knowledge in identifying COVID-19 positive and negative cases. Overall, we believe the simplicity and effectiveness of COVID-Net USPro make it a promising tool to aid the COVID-19 screening process using ultrasound images. We hope the open-source release of COVID-Net USPro helps researchers and clinical data scientists to accelerate innovations in the combat against the COVID-19 pandemic that can ultimately benefit the larger society.

Several future research directions can be explored to further improve the network. First, some additional steps in data augmentation and preparation can be taken to improve the data quality and the dataset size. In this work, ultrasound images captured with linear probes were excluded. More image augmentation and preparation techniques can be experimented with to include images captured with a linear probe and increase the data size. Second, we used simple cropping to filter out the pleura region of the images. A more procedural image segmentation step could be added to include only clinically relevant areas of the images for network construction to further improve network performance from the explainability standpoint. Lastly, we used COVIDx-US, which is a public dataset that includes data of various sources and quality. Network training could be improved by only using high-quality input ultrasound data with sufficient/specific image depth, collected, and reviewed systematically. For this purpose, a data collection protocol might be required to capture ultrasound images in a standardized manner from a set of consented participants.

Author Contributions: Conceptualization, A.E. and A.W.; methodology, A.E., J.S. and A.W.; validation, J.S., A.E., A.F. and A.W.; formal analysis, J.S.; investigation, J.S., A.E., A.F. and A.W.; resources, A.E. and A.W.; data curation, J.S. and A.E.; writing—original draft preparation, J.S. and A.E.; writing—review and editing, All authors; visualization, J.S.; supervision, A.E.; funding acquisition, A.E., P.X. and A.W. All authors have read and agreed to the published version of the manuscript.

Funding: This research was funded internally by the Pandemic Response Challenge Program of the National Research Council Canada.

Institutional Review Board Statement: Not applicable.

Informed Consent Statement: Not applicable.

Data Availability Statement: In this research, we used the COVIDx-US dataset (v1.4.). COVIDx-US is an open-access benchmark dataset of COVID-19 related ultrasound imaging data and is available to the general public at NRC COVIDx-US GitHub repository https://github.com/nrc-cnrc/COVID-US, accessed on 2 February 2023.

Conflicts of Interest: The authors declare no conflict of interest.

Sensors **2023**, 23, 2621 13 of 18

Abbreviations

SOTA

The following abbreviations are used in this manuscript:

LUSS lung ultrasound score COVID-19 coronavirus disease 2019 **POCUS** point-of-care ultrasound RT-PCR reverse transcription-polymerase chain reaction LUS lung ultrasound CXR chest X-ray CTcomputed tomography ΑI artificial intelligence CNN convolutional neural network PN prototypical network

Appendix A. Complete Network Hyperparameters

Table A1. Network Structure Hyperparameters.

state-of-the-art

Encoder	Layers with Pre-Trained Params	Layers with Trained Params	Trained Parameters (2-Way)	Trained Parameters (3-Way)	Trained Parameters (4-Way)
ResNet18L1	17	1	1026	1539	2052
ResNet18L5	13	5	8,394,754	8,395,267	8,395,780
ResNet50L1	49	1	4098	6147	8196
ResNet50L4	46	4	4,466,690	4,468,739	4,470,788

Table A2. Training Hyperparameters.

Optimizer	Adam
Reduce learning rate on plateau	
Target metric	Validation loss
Starting learning rate	0.001
Number of patience	3
Early Stopping	
Target metric	Validation loss
Number of patience	5

Appendix B. Complete Experiment Results

Table A3. Complete experiment results for ResNet18L1.

2-Way Classification								
N-Shot	Accuracy	Neg Precision	Neg Recall	COVID Precision	COVID Recall			
5	0.931	0.9473	0.924	0.9384	0.938			
10	0.962	0.9666	0.961	0.9645	0.963			
20	0.9677	0.9767	0.9597	0.9624	0.9757			
30	0.9661	0.9703	0.9627	0.9642	0.9695			
40	0.9698	0.9724	0.9677	0.9686	0.9719			
50	0.9713	0.9703	0.9732	0.9737	0.9694			
75	0.9759	0.9796	0.9723	0.9728	0.9795			
100	0.9741	0.9747	0.9738	0.9741	0.9744			
150	0.9758	0.9756	0.9763	0.9764	0.9754			
200	0.976	0.9749	0.9757	0.9754	0.975			

Sensors **2023**, 23, 2621 14 of 18

Table A3. Cont.

				3-Way Classifi	ication				
N-shot	Accuracy	Normal Precision	Normal Recall	Non-COVID Precision	Non-COVID Recall	COVID Precision	COVID Recall		
5	0.9487	0.9449	0.972	0.9783	0.937	0.947	0.937		
10	0.976	0.9672	0.989	0.9866	0.975	0.9802	0.964		
20	0.9806	0.9791	0.9905	0.9828	0.9775	0.9828	0.9737		
30	0.9832	0.9809	0.989	0.9911	0.9768	0.9794	0.9837		
40	0.982	0.979	0.9914	0.9839	0.9799	0.9847	0.9749		
50	0.9844	0.9822	0.9879	0.9892	0.9805	0.9827	0.9847		
75	0.9843	0.98	0.9915	0.9905	0.9771	0.9832	0.9844		
100	0.9863	0.9823	0.9926	0.9914	0.9807	0.9858	0.9857		
150	0.9866	0.9834	0.991	0.9906	0.9833	0.986	9854		
200	0.987	0.9841	0.9903	0.9902	0.9841	0.9881	0.9871		
				4-Way Classifi	ication				
N-shot	Accuracy	Normal Precision	Normal Recall	Other Precision	Other Recall	Non- COVID Precision	Non- COVID Recall	COVID Precision	COVID Recall
5	0.8733	0.8885	0.908	0.8361	0.798	0.9019	0.889	0.9283	0.898
10	0.8999	0.9088	0.9205	0.8653	0.831	0.9157	0.9	0.9397	0.948
20	0.9212	0.9251	0.9513	0.897	0.849	0.916	0.9268	0.9597	0.9577
30	0.927	0.9331	0.9542	0.8962	0.8657	0.9195	0.9302	0.9676	0.9578
40	0.9316	0.9374	0.9519	0.8933	0.887	0.9368	0.9263	0.9647	0.9614
50	0.9324	0.9287	0.9624	0.9103	0.8716	0.9278	0.9346	0.9669	0.9609
75	0.9377	0.9383	0.9665	0.9138	0.8855	0.9307	0.9327	0.9709	0.9663
100	0.9369	0.9382	0.9604	0.9024	0.8931	0.9353	0.9296	0.9735	0.9643
150	0.9384	0.9411	0.9655	0.9121	0.8898	0.9273	0.9335	0.9745	0.9649
200	0.941	0.9391	0.9661	0.912	0.8903	0.9348	0.941	0.9751	0.9651

 $\textbf{Table A4.} \ Complete \ experiment \ results \ for \ ResNet 18L5.$

-				2-Way Classif	ication			
N-Shot	Accuracy	Neg Precision	Neg Recall	COVID Precision	COVID Recall			
5	0.9955	0.9942	0.998	0.9983	0.993			
10	0.9947	0.9911	0.9995	0.9995	0.99			
20	0.9957	0.994	0.9977	0.9978	0.9937			
30	0.998	0.9962	0.999	0.999	0.996			
40	0.9997	0.999	0.9994	0.9994	0.9998			
50	0.9999	0.9999	0.9997	0.9997	0.9999			
75	1	1	1	1	1			
100	0.9999	1	0.9999	1	0.9999			
150	1	1	0.9999	1	0.9999			
200	1	1	1	1	1			
				3-Way Classif	ication			
N-shot	Accuracy	Normal Precision	Normal Recall	Non-COVID Precision	Non-COVID Recall	COVID Precision	COVID Recall	
5	0.9947	0.9958	0.998	0.9958	0.999	0.9992	0.992	
10	0.9903	0.9893	0.991	0.9959	0.9925	0.9888	0.9875	
20	0.9998	0.9998	0.9999	0.9998	0.9999	0.9999	0.9995	
30	0.9999	0.9997	0.9999	0.9999	0.9997	0.9999	0.9999	
40	0.9972	0.996	0.9984	0.9983	0.9991	0.9976	0.9941	
50	0.9995	0.9995	0.9999	0.999	0.9999	0.9999	0.9985	

Sensors **2023**, 23, 2621 15 of 18

Table A4. Cont.

				3-Way Classif	fication				
75	0.9988	0.9977	0.9995	0.9992	0.9991	0.9995	0.9978		
100	0.9998	0.9999	1	1	1	0.9994	0.9998		
150	1	1	1	1	1	1	1		
200	1	0.9999	1	1	0.9999	0.9999	0.9999		
				4-Way Classif	fication				
N-shot	Accuracy	Normal Precision	Normal Recall	Other Precision	Other Recall	Non- COVID Precision	Non- COVID Recall	COVID Precision	COVID Recall
5	0.9855	0.9975	0.998	0.9904	0.957	0.9643	0.99	0.998	0.997
10	0.9819	0.9932	0.996	0.9551	0.9885	0.9949	0.946	0.9923	0.997
20	0.9809	0.9796	0.9977	0.9867	0.9383	0.9623	0.9882	0.999	0.9992
30	0.9852	0.9923	0.999	0.9707	0.979	0.9815	0.963	0.999	0.9997
40	0.9881	0.999	0.999	0.999	0.954	0.9565	0.999	0.999	0.9992
50	0.983	0.9995	0.9887	0.9539	0.9882	0.9883	0.9589	0.9931	0.9962
75	0.988	0.999	0.999	0.9979	0.954	0.9565	0.998	0.999	0.999
100	0.9884	1	1	1	0.9537	0.9562	0.999	1	1
150	0.9869	1	1	0.9694	0.9791	0.9788	0.9688	1	0.9997
200	0.9891	1	1	1	0.971	0.9731	0.998	1	0.9998

 $\textbf{Table A5.} \ Complete \ experiment \ results \ for \ ResNet 50 L1.$

				2-Way Classifi	ication				
N-Shot	Accuracy	Neg Precision	Neg Recall	COVID Precision	COVID Recall				
5	0.9605	0.9784	0.947	0.9567	0.974				
10	0.974	0.9773	0.9725	0.975	0.9755				
20	0.9864	0.9854	0.9882	0.9889	0.9845				
30	0.9905	0.9919	0.9893	0.9897	0.9917				
40	0.9939	0.994	0.9939	0.994	0.9939				
50	0.9926	0.9938	0.9914	0.9916	0.9937				
75	0.9952	0.9948	0.9958	0.9959	0.9947				
100	0.9953	0.9953	0.9954	0.9954	0.9953				
150	0.9962	0.9964	0.9961	0.9961	0.9964				
200	0.9965	0.9962	0.9967	0.9967	0.9962				
				3-Way Classifi	ication				
N-shot	Accuracy	Normal Precision	Normal Recall	Non-COVID Precision	Non-COVID Recall	COVID Precision	COVID Recall		
5	0.969	0.9604	0.975	0.9847	0.971	0.9761	0.961		
10	0.9895	0.9887	0.9925	0.9933	0.9895	0.9895	0.9865		
20	0.9931	0.992	0.9935	0.9957	0.9957	0.9926	0.99		
30	0.995	0.9942	0.9962	0.996	0.9967	0.9953	0.9922		
40	0.9957	0.9959	0.9936	0.9982	0.9975	0.9933	0.9959		
50	0.9966	0.9958	0.9967	0.9981	0.9978	0.9961	0.9953		
75	0.9986	0.9968	0.9997	0.9995	0.9981	0.9996	0.9981		
100	0.9973	0.9981	0.9947	0.9993	0.9992	0.9947	0.998		
150	0.9984	0.9979	0.9974	0.9999	0.999	0.9974	0.9987		
200	0.9978	0.9981	0.998	0.9999	0.9994	0.9961	0.9971		
				4-Way Classifi	ication				
N-shot	Accuracy	Normal Precision	Normal Recall	Other Precision	Other Recall	Non- COVID Precision	Non- COVID Recall	COVID Precision	COVID Recall
5	0.932	0.9471	0.953	0.9286	0.88	0.9305	0.94	0.9622	0.955
10	0.9647	0.9789	0.98	0.9552	0.9425	0.9514	0.963	0.9869	0.9735

Sensors **2023**, 23, 2621 16 of 18

Table A5. Cont.

				4-Way Classifi	cation				
20	0.9736	0.9857	0.985	0.9651	0.9475	0.9585	0.9715	0.9901	0.9905
30	0.977	0.985	0.9892	0.9607	0.9625	0.9737	0.9642	0.9917	0.9922
40	0.9791	0.9885	0.995	0.9738	0.9546	0.9618	0.9749	0.9945	0.9917
50	0.9796	0.992	0.9934	0.978	0.948	0.9563	0.9815	0.9941	0.9955
75	0.9812	0.9935	0.9937	0.9649	0.9677	0.973	0.9675	0.9945	0.9961
100	0.9811	0.994	0.9925	0.9704	0.9601	0.9656	0.9756	0.9952	0.9962
150	0.9823	0.9946	0.9926	0.9813	0.9536	0.9591	0.9867	0.9953	0.9963
200	0.9832	0.9943	0.9928	0.9833	0.9541	0.961	0.9859	0.9951	0.9968

 Table A6. Complete experiment results for ResNet50L4.

N-shot Accuracy Neg Precision Neg Recall COVID Precision Precision COVID Recall 5 0.991 0.9919 0.992 0.9933 0.99 10 0.9955 0.999 0.991 0.999 0.9918 20 0.9982 0.9995 0.997 0.9995 0.9995 30 0.9997 0.9998 0.9995 0.9999 0.9998 40 0.9998 0.9996 0.9999 0.9999 0.9996 50 0.9982 0.9985 0.9979 0.9998 0.9998 75 0.9999 0.9999 0.9998 0.9999 0.9999 100 0.9999 0.9999 0.9999 1 0.9991 200 0.9981 0.999 0.9991 0.9999 0.9999 200 0.9999 0.9999 0.9999 0.9999 200 0.9999 0.9999 0.9999 0.9999 200 0.9999 0.9999 0.9999 0.9999	
10 0.9955 0.999 0.991 0.999 0.9918 20 0.9982 0.9995 0.997 0.9971 0.9995 30 0.9997 0.9998 0.9995 0.9995 0.9998 40 0.9998 0.9996 0.9999 0.9999 0.9996 50 0.9982 0.9985 0.9979 0.9979 0.9998 75 0.9999 0.9999 0.9998 0.9999 0.9999 100 0.9999 0.9999 0.9999 1 150 0.9981 0.999 0.9971 0.9971 0.999 200 0.999 0.999 0.999 0.999 0.999 200 0.999 0.999 0.999 0.999 0.999 N-shot Accuracy Normal Precision Precision Precision Recall Precision Recall Precision Precision Recall Precision Precision Precision Precision Precision Precision Precision Precision No.994 0.994 0.994	
20 0.9982 0.9995 0.997 0.9971 0.9995 30 0.9997 0.9998 0.9995 0.9995 0.9998 40 0.9998 0.9996 0.9999 0.9999 0.9996 50 0.9982 0.9985 0.9979 0.9985 0.9998 75 0.9999 0.9999 0.9998 0.9999 1 100 0.9999 0.9999 0.9999 1 150 0.9981 0.999 0.9991 0.9991 0.999 200 0.999 0.999 0.999 0.999 0.999 3-Way Classification N-shot Accuracy Normal Precision Precision Normal Recall Precision Recall Precision Recall Precision Recall Precision Precision Non-COVID Recall Recall Non-COVID Recall Precision Non-COVID Recall Recall Non-COVID Recall Recall Precision Non-COVID Recall Re	
30 0.9997 0.9998 0.9995 0.9995 0.9998 40 0.9998 0.9996 0.9999 0.9999 0.9996 50 0.9982 0.9985 0.9979 0.9979 0.9985 75 0.9999 0.9999 0.9998 0.9999 0.9999 100 0.9999 0.9999 0.9999 1 0.999 150 0.9981 0.999 0.9971 0.999 0.999 200 0.999 0.999 0.999 0.999 0.999 N-shot Accuracy Normal Precision Precision Normal Recall Precision Recall Precision Recall Precision Precision Precision Recall Precision No.994 0.994 0.994	
40 0.9998 0.9996 0.9999 0.9999 0.9996 50 0.9982 0.9985 0.9979 0.9979 0.9985 75 0.9999 0.9999 0.9998 0.9999 0.9999 100 0.9999 0.9999 0.9999 1 150 0.9981 0.999 0.9971 0.999 200 0.999 0.999 0.999 0.999 200 0.999 0.999 0.999 0.999 3-Way Classification N-shot Accuracy Normal Precision Precision Normal Recall Precision Recall Precision Recall Precision Precision Precision Precision Precision Normal Recall Precision Norma	
50 0.9982 0.9985 0.9979 0.9979 0.9985 75 0.9999 0.9999 0.9998 0.9999 100 0.9999 0.9999 0.9999 1 150 0.9981 0.999 0.9971 0.999 200 0.999 0.999 0.999 0.999 200 0.999 0.999 0.999 0.999 3-Way Classification N-shot Accuracy Normal Precision Precision Recall Precision Recall Precision Precision Recall Precision Recall Precision Precision Precision Precision Precision No.994 0.994	
75 0.9999 0.9998 0.9998 0.9999 100 0.9999 0.9999 0.9999 1 150 0.9981 0.999 0.9971 0.9999 200 0.999 0.999 0.999 0.999 3-Way Classification N-shot Accuracy Normal Precision Precision Normal Recall Precision Recall Precision Recall Precision Precision Recall Precision No.994 0.994 0.994	
100 0.9999 0.9999 0.9999 1 150 0.9981 0.999 0.9971 0.9991 0.999 200 0.999 0.999 0.999 0.999 3-Way Classification N-shot Accuracy Normal Precision Normal Recall Precision Recall Precision Recall Precision Precision Recall Precision Normal Recall Precision Recall Recall Precision Recall	
150 0.9981 0.999 0.9971 0.9971 0.999 200 0.999 0.999 0.999 0.999 3-Way Classification N-shot Accuracy Normal Precision Normal Recall Precision Non-COVID Recall Precision COVID Recall Precision 5 0.9953 0.9967 0.994 0.9975 0.998 0.994 0.994	
200 0.999 0.999 0.999 0.999 3-Way Classification N-shot Accuracy Normal Precision Precision Non-COVID Precision Recall Non-COVID Precision Recall COVID Precision Precision Recall COVID Recall 5 0.9953 0.9967 0.994 0.9975 0.998 0.994 0.994	
N-shot Accuracy Normal Precision Normal Recall Non-COVID Recall Precision N	
N-shot Accuracy Normal Precision Normal Recall Non-COVID Recall Precision Non-COVID Recall Precision Solution Normal Recall Precision Non-COVID Recall Precision Non-COVID Recall Precision Solution Normal Recall Precision Non-COVID Recall Precision Non-COVID Recall Precision Solution Normal Recall Precision Non-COVID Recall Precision Non-COVID Recall Precision Normal Recall Precision Normal Recall Precision Non-COVID Recall Precision Non-COVID Recall Precision Non-COVID Recall Precision Non-COVID Recall Precision Normal Recall Precision Normal Recall Precision Non-COVID Recall Precision Non-COVID Recall Precision Non-COVID Recall Non-COVID Recall Precision Non-COVID Recall Precision Non-COVID Recall Non-COVID Re	
N-Shot Accuracy Precision Normal Recall Precision COVID Recall 5 0.9953 0.9967 0.994 0.9975 0.998 0.994 0.994	
1 1111 1111 1111 1111 1111	
10 0.9958 0.9977 0.9975 0.9985 0.9915 0.9923 0.9985	
20 0.9974 0.995 0.9997 0.9998 0.998 0.9979 0.9945	
30 0.9999 0.9998 0.9998 0.9997 0.9998 0.9989 0.9979	
40 0.9937 0.9889 0.9949 0.9964 0.9985 0.9962 0.9876	
50 0.9979 0.9991 0.9955 0.9979 0.9999 0.9967 0.9981	
75 1 1 1 0.9999 1 0.9999	
100 0.9998 1 0.9994 1 1 0.9995 1	
150 1 0.9999 1 0.9999 1 1	
200 1 1 1 1 1 1 1	
4-Way Classification	
N-shot Accuracy Normal Precision Normal Recall Other Precision Other Recall Precision Recall Precision COVID	VID Recall
5 0.9777 0.9942 0.979 0.9456 0.986 0.9913 0.955 0.9967	0.991
10 0.9789 0.994 0.982 0.9696 0.966 0.9682 0.9715 0.9923	0.996
20 0.9873 0.9999 0.9998 0.9702 0.9812 0.9819 0.9682 0.9999	0.9997
30 0.9875 0.9997 0.9998 0.9926 0.958 0.9611 0.9923 0.9998	0.9999
40 0.9872 0.9995 0.9966 0.9784 0.9719 0.9737 0.9805 0.9988	0.9996
50 0.989 0.999 0.999 0.9967 0.9593 0.9616 0.9968 0.9999	0.7770
75 0.9893 0.9999 0.9989 0.9952 0.9631 0.9648 0.9953 0.9989	0.9999
100 0.9899 1 1 0.9854 0.974 0.9746 0.9854 1	
150 0.9904 1 1 0.9815 0.9802 0.9803 0.9814 1	0.9999
200 0.9908 1 1 0.9861 0.9798 0.9805 0.9824 1	0.9999

Sensors **2023**, 23, 2621 17 of 18

References

1. Ebadi, A.; Xi, P.; MacLean, A.; Florea, A.; Tremblay, S.; Kohli, S.; Wong, A. COVIDx-US: An Open-Access Benchmark Dataset of Ultrasound Imaging Data for AI-Driven COVID-19 Analytics. *Front. Biosci.-Landmark* **2022**, 27, 198. [CrossRef] [PubMed]

- 2. Cascella, M.; Rajnik, M.; Aleem, A.; Dulebohn, S.C.; Napoli, R.D. Features, Evaluation, and Treatment of Coronavirus (COVID-19). In *StatPearls*; StatPearls Publishing: Treasure Island, FL, USA, 2022.
- 3. Jacqueline, D.; Jonathan J, D.; Sarah, B.; Melissa, T.; Ada, A.; Clare, D.; Sabine, D.; Devy, E.; Yemisi, T.; Jane, C.; et al. Rapid, point-of-care antigen and molecular-based tests for diagnosis of SARS-CoV-2 infection. *Cochrane Database Syst. Rev.* **2021**, 104, 29–30. [CrossRef]
- 4. Yang, D.; Martinez, C.; Visuña, L.; Khandhar, H.; Bhatt, C.; Carretero, J. Detection and analysis of COVID-19 in medical images using Deep Learning Techniques. *Sci. Rep.* **2021**, *11*, 19638. . [CrossRef] [PubMed]
- 5. Wang, L.; Lin, Z.Q.; Wong, A. COVID-Net: A tailored deep convolutional neural network design for detection of COVID-19 cases from chest X-ray images. *Sci. Rep.* **2020**, *10*, 19549. [CrossRef] [PubMed]
- 6. Gunraj, H.; Sabri, A.; Koff, D.; Wong, A. COVID-Net CT-2: Enhanced Deep Neural Networks for Detection of COVID-19 From Chest CT Images Through Bigger, More Diverse Learning. *Front. Med.* **2022**, *8*, 729287. [CrossRef] [PubMed]
- 7. MacLean, A.; Abbasi, S.; Ebadi, A.; Zhao, A.; Pavlova, M.; Gunraj, H.; Xi, P.; Kohli, S.; Wong, A. COVID-Net US: A Tailored, Highly Efficient, Self-Attention Deep Convolutional Neural Network Design for Detection of COVID-19 Patient Cases from Point-of-care Ultrasound Imaging. In Proceedings of the FAIR-MICCAI'21, Strasbourg, France, 27 September–1 October 2021.
- 8. Kasban, H. A Comparative Study of Medical Imaging Techniques. Int. J. Inf. Sci. Intell. Syst. 2015, 4, 37–58.
- Gunraj, H.; Wang, L.; Wong, A. COVIDNet-CT: A Tailored Deep Convolutional Neural Network Design for Detection of COVID-19 Cases From Chest CT Images. Front. Med. 2020, 7, 1025. [CrossRef] [PubMed]
- 10. Snell, J.; Swersky, K.; Zemel, R. Prototypical Networks for Few-shot Learning. In Proceedings of the Advances in Neural Information Processing Systems, Long Beach, CA, USA, 4–9 December 2017; Guyon, I., Luxburg, U.V., Bengio, S., Wallach, H., Fergus, R., Vishwanathan, S., Garnett, R., Eds.; Curran Associates, Inc.: Red Hook, NY, USA, 2017; Volume 30.
- 11. Li, G.; Jampani, V.; Sevilla-Lara, L.; Sun, D.; Kim, J.; Kim, J. Adaptive Prototype Learning and Allocation for Few-Shot Segmentation. In Proceedings of the IEEE/CVF Conference on Computer Vision and Pattern Recognition (CVPR), Nashville, TN, USA, 20–25 June 2021; pp. 8334–8343.
- 12. Sun, S.; Sun, Q.; Zhou, K.; Lv, T. Hierarchical Attention Prototypical Networks for Few-Shot Text Classification. In Proceedings of the 2019 Conference on Empirical Methods in Natural Language Processing and the 9th International Joint Conference on Natural Language Processing (EMNLP-IJCNLP), Hong Kong, China, 3–7 November 2019; Association for Computational Linguistics: Hong Kong, China, 2019; pp. 476–485. [CrossRef]
- 13. Deuschel, J.; Firmbach, D.; Geppert, C.I.; Eckstein, M.; Hartmann, A.; Bruns, V.; Kuritcyn, P.; Dexl, J.; Hartmann, D.; Perrin, D.; et al. Multi-Prototype Few-Shot Learning in Histopathology. In Proceedings of the IEEE/CVF International Conference on Computer Vision (ICCV) Workshops, Montreal, BC, Canada, 11–17 October 2021; pp. 620–628.
- 14. Selvaraju, R.R.; Cogswell, M.; Das, A.; Vedantam, R.; Parikh, D.; Batra, D. Grad-CAM: Visual Explanations from Deep Networks via Gradient-Based Localization. In Proceedings of the 2017 IEEE International Conference on Computer Vision (ICCV), Venice, Italy, 22–29 October 2017; pp. 618–626. [CrossRef]
- 15. Lin, Z.Q.; Shafiee, M.J.; Bochkarev, S.; Jules, M.S.; Wang, X.; Wong, A. Do Explanations Reflect Decisions? A Machine-centric Strategy to Quantify the Performance of Explainability Algorithms. *arXiv* **2019**, arXiv:1910.07387.
- 16. Wong, A. COVID-Net Open Initiative. Available online: https://alexswong.github.io/COVID-Net/ (accessed on 2 February 2023).
- 17. Aboutalebi, H.; Pavlova, M.; Gunraj, H.; Shafiee, M.J.; Sabri, A.; Alaref, A.; Wong, A. MEDUSA: Multi-scale Encoder-Decoder Self-Attention Deep Neural Network Architecture for Medical Image Analysis. Front. Med. 2022, 8, 2891. [CrossRef] [PubMed]
- 18. Ozturk, T.; Talo, M.; Yildirim, E.A.; Baloglu, U.B.; Yildirim, O.; Rajendra Acharya, U. Automated detection of COVID-19 cases using deep neural networks with X-ray images. *Comput. Biol. Med.* **2020**, *121*, 103792. [CrossRef] [PubMed]
- 19. Afshar, P.; Heidarian, S.; Naderkhani, F.; Oikonomou, A.; Plataniotis, K.N.; Mohammadi, A. COVID-CAPS: A capsule network-based framework for identification of COVID-19 cases from X-ray images. *Pattern Recognit. Lett.* **2020**, *138*, 638–643. [CrossRef] [PubMed]
- 20. Thakur, S.; Kumar, A. X-ray and CT-scan-based automated detection and classification of covid-19 using convolutional neural networks (CNN). *Biomed. Signal Process. Control* **2021**, *69*, 102920. [CrossRef]
- 21. Diaz-Escobar, J.; Ordóñez-Guillén, N.E.; Villarreal-Reyes, S.; Galaviz-Mosqueda, A.; Kober, V.; Rivera-Rodriguez, R.; Lozano Rizk, J.E. Deep-learning based detection of COVID-19 using lung ultrasound imagery. *PLoS ONE* **2021**, *16*, e0255886. [CrossRef] [PubMed]
- 22. Simonyan, K.; Zisserman, A. Very Deep Convolutional Networks for Large-Scale Image Recognition. In Proceedings of the International Conference on Learning Representations, San Diego, CA, USA, 7–9 May 2015.
- 23. Szegedy, C.; Vanhoucke, V.; Ioffe, S.; Shlens, J.; Wojna, Z. Rethinking the Inception Architecture for Computer Vision. In Proceedings of the IEEE Conference on Computer Vision and Pattern Recognition, Las Vegas, NV, USA, 26 June–1 July 2016.
- 24. He, K.; Zhang, X.; Ren, S.; Sun, J. Deep Residual Learning for Image Recognition. In Proceedings of the IEEE Conference on Computer Vision and Pattern Recognition, Las Vegas, NV, USA, 26 June–1 July 2016.

Sensors **2023**, 23, 2621 18 of 18

 Shorfuzzaman, M.; Hossain, M.S. MetaCOVID: A Siamese neural network framework with contrastive loss for n-shot diagnosis of COVID-19 patients. *Pattern Recognit.* 2021, 113, 107700. [CrossRef] [PubMed]

- 26. Ebadi, A.; Azimi, H.; Xi, P.; Tremblay, S.; Wong, A. COVID-Net FewSE: An Open-Source Deep Siamese Convolutional Network Model for Few-Shot Detection of COVID-19 Infection from X-ray Images. *J. Comput. Vis. Imaging Syst.* **2021**, *7*, 16–18.
- 27. Karnes, M.; Perera, S.; Adhikari, S.; Yilmaz, A. Adaptive Few-Shot Learning PoC Ultrasound COVID-19 Diagnostic System. In Proceedings of the 2021 IEEE Biomedical Circuits and Systems Conference (BioCAS), Virtual, 7–9 October 2021. [CrossRef]
- 28. Lee, W.; Roh, Y. Ultrasonic transducers for medical diagnostic imaging. Biomed. Eng. Lett. 2017, 7, 91–97. [CrossRef]
- 29. Hussain, Z.; Gimenez, F.; Yi, D.; Rubin, D. Differential data augmentation techniques for medical imaging classification tasks. In Proceedings of the AMIA Annual Symposium Proceedings 2018, San Francisco, CA, USA, 3–7 November 2018.
- 30. Amann, J.; Blasimme, A.; Vayena, E.; Frey, D.; Madai, V.I. Explainability for Artificial Intelligence in healthcare: A multidisciplinary perspective. *BMC Med. Inform. Decis. Mak.* **2020**, 20, 310. [CrossRef] [PubMed]
- 31. Huang, G.; Liu, Z.; Pleiss, G.; Van Der Maaten, L.; Weinberger, K. Convolutional Networks with Dense Connectivity. *IEEE Trans. Pattern Anal. Mach. Intell.* **2019**, 44, 8704–8716. [CrossRef] [PubMed]
- 32. Deng, J.; Dong, W.; Socher, R.; Li, L.J.; Li, K.; Fei-Fei, L. Imagenet: A large-scale hierarchical image database. In Proceedings of the 2009 IEEE Conference on Computer Vision and Pattern Recognition, Miami, FL, USA, 20–25 June 2009; IEEE: Piscataway, NJ, USA, 2009; pp. 248–255.
- 33. Lichtenstein, D.A.; Mezière, G.A.; Lagoueyte, J.F.; Biderman, P.; Goldstein, I.; Gepner, A. A-lines and B-lines: Lung ultrasound as a bedside tool for predicting pulmonary artery occlusion pressure in the critically ill. *Chest* **2009**, *136*, 1014–1020. [CrossRef] [PubMed]

Disclaimer/Publisher's Note: The statements, opinions and data contained in all publications are solely those of the individual author(s) and contributor(s) and not of MDPI and/or the editor(s). MDPI and/or the editor(s) disclaim responsibility for any injury to people or property resulting from any ideas, methods, instructions or products referred to in the content.