

## NRC Publications Archive Archives des publications du CNRC

### **Binding symmetry and surface flexibility mediate antibody self-association**

Schrag, Joseph D.; Picard, Marie-Ève; Gaudreault, Francis; Gagnon, Louis-Patrick; Baardsnes, Jason; Manenda, Mahder S.; Sheff, Joey; Deprez, Christophe; Baptista, Cassio; Hogues, Hervé; Kelly, John F.; Purisima, Enrico O.; Shi, Rong; Sulea, Traian

This publication could be one of several versions: author's original, accepted manuscript or the publisher's version. / La version de cette publication peut être l'une des suivantes : la version prépublication de l'auteur, la version acceptée du manuscrit ou la version de l'éditeur.

For the publisher's version, please access the DOI link below. / Pour consulter la version de l'éditeur, utilisez le lien DOI ci-dessous.

#### **Publisher's version / Version de l'éditeur:**

<https://doi.org/10.1080/19420862.2019.1632114>

*mAbs*, 11, 7, pp. 1300-1318, 2019-07-19

#### **NRC Publications Archive Record / Notice des Archives des publications du CNRC :**

<https://nrc-publications.canada.ca/eng/view/object/?id=0c949240-61b0-4677-8818-225315b4f12a>

<https://publications-cnrc.canada.ca/fra/voir/objet/?id=0c949240-61b0-4677-8818-225315b4f12a>

Access and use of this website and the material on it are subject to the Terms and Conditions set forth at

<https://nrc-publications.canada.ca/eng/copyright>

READ THESE TERMS AND CONDITIONS CAREFULLY BEFORE USING THIS WEBSITE.

L'accès à ce site Web et l'utilisation de son contenu sont assujettis aux conditions présentées dans le site


<https://publications-cnrc.canada.ca/fra/droits>

LISEZ CES CONDITIONS ATTENTIVEMENT AVANT D'UTILISER CE SITE WEB.

**Questions?** Contact the NRC Publications Archive team at

PublicationsArchive-ArchivesPublications@nrc-cnrc.gc.ca. If you wish to email the authors directly, please see the first page of the publication for their contact information.

**Vous avez des questions?** Nous pouvons vous aider. Pour communiquer directement avec un auteur, consultez la première page de la revue dans laquelle son article a été publié afin de trouver ses coordonnées. Si vous n'arrivez pas à les repérer, communiquez avec nous à PublicationsArchive-ArchivesPublications@nrc-cnrc.gc.ca.



## Binding symmetry and surface flexibility mediate antibody self-association

Joseph D. Schrag, Marie-Ève Picard, Francis Gaudreault, Louis-Patrick Gagnon, Jason Baardsnes, Mahder S. Manenda, Joey Sheff, Christophe Deprez, Cassio Baptista, Hervé Hogues, John F. Kelly, Enrico O. Purisima, Rong Shi & Traian Sulea

To cite this article: Joseph D. Schrag, Marie-Ève Picard, Francis Gaudreault, Louis-Patrick Gagnon, Jason Baardsnes, Mahder S. Manenda, Joey Sheff, Christophe Deprez, Cassio Baptista, Hervé Hogues, John F. Kelly, Enrico O. Purisima, Rong Shi & Traian Sulea (2019) Binding symmetry and surface flexibility mediate antibody self-association, *mAbs*, 11:7, 1300-1318, DOI: [10.1080/19420862.2019.1632114](https://doi.org/10.1080/19420862.2019.1632114)

To link to this article: <https://doi.org/10.1080/19420862.2019.1632114>



© 2019 The Author(s). Published with license by Taylor & Francis Group, LLC.



View supplementary material [↗](#)



Published online: 18 Jul 2019.



Submit your article to this journal [↗](#)



Article views: 1878



View related articles [↗](#)












View Crossmark data [↗](#)



Citing articles: 1 View citing articles [↗](#)

## Binding symmetry and surface flexibility mediate antibody self-association

Joseph D. Schrag<sup>\*a</sup>, Marie-Ève Picard <sup>\*b</sup>, Francis Gaudreault <sup>\*a</sup>, Louis-Patrick Gagnon <sup>a</sup>, Jason Baardsnes <sup>a</sup>, Mahder S. Manenda <sup>b</sup>, Joey Sheff <sup>c</sup>, Christophe Deprez <sup>a</sup>, Cassio Baptista<sup>a</sup>, Hervé Hogues<sup>a</sup>, John F. Kelly<sup>c</sup>, Enrico O. Purisima<sup>a</sup>, Rong Shi <sup>b</sup>, and Traian Sulea <sup>a</sup>

<sup>a</sup>Human Health Therapeutics Research Centre, National Research Council Canada, Montreal, QC H4P 2R2, Canada; <sup>b</sup>Département de Biochimie, de Microbiologie et de Bio-informatique, PROTEO, and Institut de Biologie Intégrative et des Systèmes (IBIS), Université Laval, Pavillon Charles-Eugène-Marchand, Québec City, QC G1V 0A6, Canada; <sup>c</sup>Human Health Therapeutics Research Centre, National Research Council Canada, Ottawa, ON K1A 0R6, Canada

### ABSTRACT

Solution stability is an important factor in the optimization of engineered biotherapeutic candidates such as monoclonal antibodies because of its possible effects on manufacturability, pharmacology, efficacy and safety. A detailed atomic understanding of the mechanisms governing self-association of natively folded protein monomers is required to devise predictive tools to guide screening and re-engineering along the drug development pipeline. We investigated pairs of affinity-matured full-size antibodies and observed drastically different propensities to aggregate from variants differing by a single amino-acid. Biophysical testing showed that antigen-binding fragments (Fabs) from the aggregating antibodies also reversibly associated with equilibrium dissociation constants in the low-micromolar range. Crystal structures (PDB accession codes 6MXR, 6MXS, 6MY4, 6MY5) and bottom-up hydrogen-exchange mass spectrometry revealed that Fab self-association occurs in a symmetric mode that involves the antigen complementarity-determining regions. Subtle local conformational changes incurred upon point mutation of monomeric variants foster formation of complementary polar interactions and hydrophobic contacts to generate a dimeric Fab interface. Testing of popular *in silico* tools generally indicated low reliabilities for predicting the aggregation propensities observed. A structure-aggregation data set is provided here in order to stimulate further improvements of *in silico* tools for prediction of native aggregation. Incorporation of intermolecular docking, conformational flexibility, and short-range packing interactions may all be necessary features of the ideal algorithm.

### ARTICLE HISTORY

Received 18 March 2019  
Revised 6 June 2019  
Accepted 11 June 2019

### KEYWORDS

Aggregation; native folding; single point mutation; structure-aggregation relationship; prediction method



## Introduction

Protein instability in solution is a long-standing problem in the biopharmaceutical industry. It affects many critical aspects of the biologic drug discovery and development process, including binding affinity, biological activity, protein expression, manufacturability, storage, handling, delivery, efficacy, and safety.<sup>1–7</sup> Controlling and mitigating aggregation problems by molecular re-engineering in such a multi-factorial landscape is a challenging task.<sup>8,9</sup> The challenge is further exacerbated by a relatively limited understanding of aggregation mechanisms at the molecular level, which impacts our ability to accurately predict aggregation propensities upon molecular redesign or modification.<sup>10–12</sup>

Protein aggregation is commonly perceived as resulting from non-specific interactions between hydrophobic core regions exposed by unfolding or misfolding and is related to the conformational or folding stability.<sup>10,13</sup> Another critical, albeit somewhat under-appreciated<sup>14</sup> and less well understood,<sup>9</sup> pathway for protein aggregation is self-association of natively folded protein monomers.<sup>4</sup> These interactions typically lead to


reversible homo-oligomers that may nucleate irreversible aggregates in dilute solution,<sup>10</sup> and contributes to viscosity at high concentrations required for sub-cutaneous administration.<sup>3,15</sup> In dilute solution, native self-association is thought to occur mainly due to exposure of hydrophobic patches at the surface of the folded monomers, with some contributions from domain-level charge distributions,<sup>11,16,17</sup> whereas at high concentration it is mainly driven by electrostatic forces between the charge distributions of the folded structures with secondary contributions from hydrophobic interactions.<sup>18–21</sup>

Much of our understanding of protein-protein interactions and self-assembly at the molecular level has been derived from analysis of crystal structures and crystal packing, and more recently from cryo-electron microscopy (cryo-EM) studies.<sup>22,23</sup> Substantial efforts have been dedicated to characterizing the self-assembling interfaces deduced from crystal structures in order to discern specific homodimers existing in solution from packing interactions occurring at high concentration during crystallization.<sup>24–29</sup> This included calculations of interfacial properties such as size, shape, spatial distribution, chemical composition, hydrogen

**CONTACT** Traian Sulea  [traian.sulea@nrc-cnrc.gc.ca](mailto:traian.sulea@nrc-cnrc.gc.ca)  Human Health Therapeutics Research Centre, National Research Council Canada, 6100 Royalmount Ave., Montreal, QC H4P 2R2, Canada

Color versions of one or more of the figures in the article can be found online at [www.tandfonline.com/kmab](http://www.tandfonline.com/kmab).

\*These authors contributed equally to this work.

 Supplemental data for this article can be accessed on the [publisher's website](#).

© 2019 The Author(s). Published with license by Taylor & Francis Group, LLC.

This is an Open Access article distributed under the terms of the Creative Commons Attribution-NonCommercial-NoDerivatives License (<http://creativecommons.org/licenses/by-nc-nd/4.0/>), which permits non-commercial re-use, distribution, and reproduction in any medium, provided the original work is properly cited, and is not altered, transformed, or built upon in any way.

bonding, shape complementarity and interfacial water, with certain property combinations especially helpful in detecting true self-association. In general, larger and more hydrophobic interfaces with good shape complementarity and fewer internal cavities filled with water are indicative of specific oligomeric interactions, which can be accompanied by local induced-fit conformational changes. The development of the core-rim spatial architecture of protein-protein interfaces also indicated a “dry” hydrophobic interface-core surrounded by a “wet” hydrophilic interface-rim in the case of specific associations, whereas no core can be distinguished for crystal contacts due to small buried area and residue composition resembling typical protein surfaces.<sup>30–32</sup> The more objective definition of the newer core-support-rim interface architecture allows calculation of interface propensity scores that can be used to predict the self-association effects of mutations within a given surface patch.<sup>33</sup>

The fundamental role of symmetric oligomerization in biochemistry was recognized long ago by the Monod-Wyman-Changeaux model of allosteric transitions of proteins made up of identical subunits.<sup>34</sup> It became clear that the difference between monomeric and oligomeric structures is small and can have large effects on structure and function. For example, a single mutation in hemoglobin leads to filament formation and sickle cell.<sup>35</sup> The importance of symmetry in oligomerization was underscored by studies in which symmetry was used to design oligomerization interfaces and by random mutations that led to assembly into fibers.<sup>36</sup> Perhaps one of the most compelling molecular views of reversible self-interactions between natively folded proteins is supramolecular symmetric self-assembly into fibers gleaned from single-particle cryo-EM at 7–8 Å resolution.<sup>14</sup> One fundamental observation stemming from this study is that the naturally evolved tendency of proteins towards homo-oligomerization is deeply rooted in symmetry, which amplifies energetic effects,<sup>37</sup> triggering formation or destruction of large self-assemblies dependent on a single mutation.

In two cases during the development of therapeutic antibodies, X-ray crystallography has been employed to associate aggregation or high viscosity with symmetric homo-dimerization of the antigen-binding fragments (Fabs) *via* the antigen complementarity-determining region (CDR),<sup>38,39</sup> and in a third case it has been associated with symmetric homo-tetramerization of the Fabs, also implicating CDR residues.<sup>40–42</sup> In these cases, a few mutations at the crystallographically observed self-association interface sufficed to eliminate aggregation, suggesting that the contacts observed in the crystal lattice were not merely due to crystal packing, but may reflect true self-assembly modes in solution. Some success in preventing aggregation *via* mutagenesis targeting aggregation hot spots mapped within the CDR has also been reported, although the implicated self-association interfaces were not directly elucidated at the atomic level.<sup>9</sup>


We recently conducted a structure-based affinity maturation campaign in which antigen binding and not self-association represented the pressure along an “evolutionary” pathway consisting of only a few dozens of variants.<sup>43</sup> Affinity-matured IgG variants differing by a single point mutation were observed to differ substantially in their propensity to aggregate. In this study, we describe biophysical and structural characterization of these variants and their Fabs and investigate the structural determinants of the self-association. We also characterized the aggregation propensities of

the other mutants generated along the original affinity maturation campaign, and built a Structure-native-Aggregation-Relationship (SnAggR) dataset around the crystal structures determined for self-associating Fab variants. This dataset deepens our mechanistic understanding of self-association between natively folded monomers, and could spur improvements in *in silico* structure-based prediction methods for native-folded aggregation.<sup>11,16,17</sup>

## Results

### Single substitutions lead to distinct aggregation profiles of full-length antibodies

The parental anti-vascular endothelial growth factor (VEGF)-A Fab bH1,<sup>44,45</sup> and eight variants affinity-matured *via* the Assisted Design of Antibody and Protein Therapeutics (ADAPT) platform,<sup>43</sup> were produced as full-length human IgG1 antibodies and underwent a developability assessment in terms of aggregation propensity. The variants included four triple heavy chain mutants (listed in Figure 1) and four variants containing the same triple heavy chain mutations plus an additional light chain mutation. The parental and four affinity-matured variants behave mainly as monomeric antibodies. By size exclusion chromatography (SEC) we observed a large peak accounting for ~90% of the absorbance with an apparent molecular weight of ~150 kDa, as expected for monomeric protein (Figure S1). Small peaks accounting for the remainder of the absorbance have apparent molecular weights consistent with small oligomers, most of which are expected to be dimers and trimers. Dynamic light scattering (DLS) analysis of the major SEC peak indicates particles with a hydrodynamic radius of ~5.5 nm and a self-diffusion coefficient,  $D_0$ , of  $\sim 4.5 \times 10^{-7}$  cm<sup>2</sup>/s (Figure 2(a,b), Table S1) after 2 weeks of storage at 40°C, both values being consistent with monomeric IgG.<sup>46,47</sup> Similar diffusion data were obtained at 4°C (Table S1). Sedimentation velocity analytical ultracentrifugation (AUC) data showed a  $c(s)$  distribution with >90% of the population sedimenting as 6.5 S particles (Figure 2(c)), consistent with monomeric IgG. Analysis of the minor SEC peak showed a predominance of ~9 S particles with a hydrodynamic radius of ~8 nm and a diffusion coefficient of  $\sim 3 \times 10^{-7}$  cm<sup>2</sup>/s, all consistent with dimeric IgG. In contrast, the other four affinity matured-variants eluted from SEC columns as broad peaks with unusual elution times indicative of stationary phase interactions, had a broad distribution of hydrodynamic radii greater than 10 nm (Figure 2(a)), diffusion coefficients of  $1\text{--}2 \times 10^{-7}$  cm<sup>2</sup>/s (Figure 2(b), Table S1), and showed a broad  $c(s)$  distribution with a predominance of 12–15 S particles (Figure 2(c)), all suggesting aggregation into larger multimers, perhaps ranging from trimers to pentamers. These variants remained oligomers at concentrations ranging from 1–20 mg/mL. The melting temperatures ( $T_m$ ) of all variants as measured by differential scanning calorimetry (DSC) were found to be similar to that of the parental full-length antibody with a major transition at ~72–76°C and a secondary transition around 83°C (Figure 2(d)), suggesting that self-interactions of the aggregating variants are not a result of major conformational changes that alter thermal stability. The difference between aggregating affinity-matured triple mutant monoclonal antibodies (mAbs) and those that did not aggregate was a single amino



**Figure 1.** Affinity-matured antibody variants with distinct aggregation behaviors due to single mutation. Location of key positions for affinity maturation on the crystal structure (PDB code 3BDY) of bH1-Fab (heavy chain shown in dark blue, light chain shown in light blue) bound to the VEGF antigen (gray). Affinity maturation data taken from ref. 43.

acid substitution, Met *versus* Asp at heavy-chain position <sup>H</sup>99 (Chothia numbering adopted throughout this work) (Figure 1).

Negative stain transmission electron microscopy (TEM) analysis confirmed these findings. The major SEC peak of the antibody variant having the G<sup>H</sup>99D mutation contains monodisperse antibodies with the expected “Y” shape (Figure 2(e), left image). By contrast, the antibody variant having the G<sup>H</sup>99M mutation shows an ensemble of oligomeric structures of varying sizes and geometries, including open-ended as well as circular shapes with discernable “Y”-shaped building blocks (Figure 2(e), right image).


### Fabs of aggregating antibodies form dimers

In order to shed light into the aggregation mechanism of the full-length IgG variants, we analyzed the self-association of the corresponding Fabs of the affinity-matured variants. At a concentration of 1 mg/mL, sedimentation velocity data indicated over 97% monomer in the cases of the parental bH1 Fab and the affinity-matured triple mutants bearing Asp at position <sup>H</sup>99, whereas about 90% dimer formation was observed for the affinity-matured variants with Met at this position (Figure 3(a)). This strongly suggests that the full-length antibody variants self-associate *via* Fab-Fab intermolecular interactions. The equilibrium dissociation constant ( $K_D$ ) calculated from the sedimentation velocity data using a monomer-dimer fast equilibrium model gave values of 10–50  $\mu\text{M}$ . By isothermal titration calorimetry (ITC) the  $K_D$  values for the Fab variants forming dimers were determined to be below 10  $\mu\text{M}$  with a large exothermic enthalpy contribution, whereas no binding was detected for the monomeric Fab variants (Figure 3(b)). DLS experiments further corroborated the AUC and ITC data, with lower intrinsic diffusion coefficients ( $D_0$ ) for the preferentially dimeric variants relative to the predominantly monomeric ones after two weeks of storage at 40°C (Figure 3(c,d), Table S2). Similar DLS data were obtained at 4°C

(Table S2), especially in terms of self-diffusion coefficients,  $D_0$ . Diffusion interaction parameter,  $k_D$ , values were generally more negative after two-weeks storage at 40°C than with storage at 4°C, reflecting an increased tendency for self-interactions due to the stress imposed by storage at higher temperature. The thermal stabilities of these Fab variants are similar to those of the parental bH1-Fab, with  $T_m$  values measured by DSC around 76°C (Figure 3(e)), which further suggests self-assembly between native-like folded monomers and not *via* an unfolding pathway.<sup>4</sup>

### Fab dimers display symmetric geometries

In order to gain an atomic-level understanding of Fab dimerization leading to the aggregation of full-length antibodies, we determined crystal structures of the four dimeric Fab variants. In each case the asymmetric unit of the crystal is a dimer. As shown in Figure 4 for the Fab variant W<sup>H</sup>33-F<sup>H</sup>98-M<sup>H</sup>99, dimerization occurs *via* a non-crystallographic 2-fold symmetry axis mediating head-to-head interactions of the CDRs. The total solvent-accessible surface area (SASA) of two Fab molecules buried upon their association is approximately 2600  $\text{\AA}^2$ , with a per-monomer interfacial area of 1294  $\text{\AA}^2$ . This is well within the range observed for protein-protein interfaces,<sup>28</sup> and is larger than typical antigen-antibody interfaces,<sup>26</sup> including the parental bH1 Fab-VEGF complex, which has a total buried SASA of 1542  $\text{\AA}^2$ . Distinguishing dimer interfaces from crystal contacts based on buried surface areas is imperfect, but the interfacial areas observed here are more typical of homodimer interactions than crystal contacts.<sup>24,25,27,29</sup> In addition to buried interface area, we calculated a dozen interface properties proposed to help discern true homodimeric interface in solution from crystal packing (Table S3). Taken together, these data confirm that the observed symmetric dimer interface does not fit the crystal packing characteristics. In particular, two powerful discriminating property combinations consisting of either: 1) buried



**Figure 2.** Aggregation behavior of parental bH1 and affinity-matured human IgG1/k full-size antibody variants. (a) DLS intensity distributions. (b) Concentration dependence of the diffusion coefficient as measured by DLS. (c) Sedimentation velocity  $c(s)$  distributions normalized to equal total areas. (d) DSC thermograms. In A-D solid lines or triangles represent the main peaks from SEC of parental bH1 (black lines and symbols) or variants containing D<sup>H</sup>99. The minor SEC peaks of parental bH1 and variants containing D<sup>H</sup>99 are represented with broken lines (---) or circles. Dashed lines (---) or diamonds represent variants with M<sup>H</sup>99. (e) Negative stain transmission electron microscopy (TEM) images for a D<sup>H</sup>99 containing variant (left) and an M<sup>H</sup>99 (right) affinity-matured variant. Sample image at 280k magnification with a 50-nm scale bar. Inset shows 2D class average of independent particles.

nonpolar interface area ( $f_{np}B$ ) of 1294 Å<sup>2</sup>, fraction of buried atoms ( $f_{bu}$ ) of 41% and residue propensity (RP) score of 3.7;<sup>24</sup> or 2) buried interface area ( $B$ ) of 2588 Å<sup>2</sup> and gap volume (GV) index of 0.75<sup>25</sup> assign this interface to the soluble homodimeric class. Crystal structures determined for W<sup>H</sup>33-M<sup>H</sup>98-M<sup>H</sup>99 and


two other dimeric Fab variants indicate the same mode of assembly (Table 1 and Figure S2) with similar interface metrics (Table S3). Careful analysis of the electron density maps showed no evidence of methionine oxidation at position <sup>H</sup>99, further supported by comparison with electron density maps of high-

resolution structures ( $< 1.4 \text{ \AA}$ ) from the Protein Data Bank (PDB) containing methionine sulfoxides. Attempts to crystallize a monomeric Fab variant, the parental bH1, were unsuccessful, possibly due to the absence of significant Fab self-interactions even at the high concentration used for crystallization.

### HX-MS supports symmetric self-assembly configuration

Hydrogen exchange – mass spectrometry (HX-MS) was used to confirm the dimerization mode in solution. HX-MS has previously been used to map conformational behavior of self-associating mAb variants.<sup>48,49</sup> The high propensity for dimerization at workable experimental concentrations of the  $W^{H33}\text{-F}^{H98}\text{-M}^{H99}$  variant and strictly monomeric behavior of the  $W^{H33}\text{-F}^{H98}\text{-D}^{H99}$  variant prevented the direct comparison of distinct molecular assemblies for each variant. Therefore, the summed deuteration

( $D_s$ ) over 5 timepoints spanning 0.05 to 60 min was compared between the monomeric ( $W^{H33}\text{-F}^{H98}\text{-D}^{H99}$ ) and the dimeric ( $W^{H33}\text{-F}^{H98}\text{-M}^{H99}$ ) Fab variants (Figure S3). The full set of kinetic plots are shown in Figure S4. A total of 57 heavy-chain and 57 light-chain peptides were analyzed, resulting in 80 and 99% total sequence coverage, respectively. Peptides overlapping position  $H^{99}$  were excluded due to inherent differences in the intrinsic exchange rates of Asp and Met mutants, accounting for the lower overall coverage of the heavy chain. A reduction in deuteration was observed in heavy-chain peptides spanning residues  $N^{H28}\text{-F}^{H67}$  (CDR-H1/2) and  $Y^{H100a}\text{-L}^{H110}$  (CDR-H3), and light-chain peptides  $C^{L23}\text{-G}^{L31}$  (CDR-L1),  $L^{L47}\text{-L}^{L73}$  (CDR-L2), and  $Y^{L87}\text{-F}^{L115}$  (CDR-L3). Significant differences in summed deuteration ( $\Delta D_s$ ) were mostly confined to the CDR regions and the proposed dimerization interface (Figure S3C). The lack of significant differences in non-interfacial regions suggests that the global




**Figure 3.** Self-association behavior of parental bH1 and affinity-matured Fab variants. (a) Sedimentation velocity  $c(s)$  distributions normalized to equal total areas. (b) ITC plots showing the integrated heats of select Fab variants titrated into buffer and curve-fit to the dimer-dissociation model. (c) DLS intensity distributions. (d) Concentration dependence of diffusion coefficients as measured by DLS. (e) DSC thermograms. Solid lines or triangles represent parental bH1 (black lines or symbols) or variants containing  $D^{H99}$ . Dashed lines or diamonds represent variants with  $M^{H99}$ .

conformational impact of the mutation is minor. Therefore, observed mass shifts can be attributed to local conformational divergence between the two Fab variants and increased stabilization as a result of dimerization. Our HX results suggest that the symmetric geometry of self-assembly of the Fab variants persists in solution and is not an artifact of crystal packing.

### Self-assembly engages short-range complementary interactions

The crystal structure of the dimeric bH1-Fab variant  $W^H33-F^H98-M^H99$  (for primary sequence information see Figure S5) solved at 1.95-Å resolution provides a detailed view of the self-assembly interface. Symmetrical engagement between the 6 CDR loops from the two Fab monomers leads to numerous polar and non-polar contacts involving predominantly side-chain atoms. No less than 12 direct hydrogen (H)-bonds can be established between



**Figure 4.** Crystal structure of the dimeric bH1-Fab variant  $W^H33-F^H98-M^H99$ , representative of other dimeric Fab structures determined in this study. (a) The overall geometry is shown on top from two viewing directions. The interface region is expanded in the bottom panels. Interacting Fab molecules are colored in shades of red and blue, with the heavy chains in darker tones and the light chains in paler tones. In the zoomed-in panels, residues establishing intermolecular contacts are rendered as sticks, H-bonds are shown as yellow dashes, and interfacial water molecules and chloride ions are shown as red and green spheres, respectively. Stereoviews detailing dimer interface interactions around key positions, including: (b) inter-Fab contacts implicating CDR-H1 residue  $W^H33$  and CDR-H3 residues  $F^H98$  and  $M^H99$ , and (c) inter-Fab contacts established by CDR-L1 residue  $R^L30b$  present in the related variant  $W^H33-F^H98-M^H99-R^L30b$ . Residues corresponding to the positions probed in the SnAggr dataset from Table 2 are highlighted by translucent sphere models. Color coding of protein chains is as in (a).

**Table 1.** X-ray data collection and refinement statistics.


Structure	W <sup>H</sup> 33-M <sup>H</sup> 98-M <sup>H</sup> 99 bH1 mutant	W <sup>H</sup> 33-F <sup>H</sup> 98-M <sup>H</sup> 99 bH1 mutant	W <sup>H</sup> 33-M <sup>H</sup> 98-M <sup>H</sup> 99-R <sup>L</sup> 30b bH1 mutant	W <sup>H</sup> 33-F <sup>H</sup> 98-M <sup>H</sup> 99-R <sup>L</sup> 30b bH1 mutant
Space group	P21	P21	P21	P21
<i>a</i> , <i>b</i> , <i>c</i> (Å)	65.09, 95.47, 109.73	65.68, 95.96, 110.11	81.87, 65.86, 93.13	81.60, 65.53, 92.56
$\alpha$ , $\beta$ , $\gamma$ (°)	90, 103.01, 90	90, 103.68, 90	90, 99.93, 90	90, 100.73, 90
Wavelength (Å)	0.97949	0.97949	0.97949	0.97949
Resolution (Å)	50.0–2.04 (2.08–2.04)	50.0–1.95 (1.98–1.95)	50.0 – 1.69 (1.72–1.69)	50–1.73 (1.76–1.73)
Observed <i>hkl</i>	281,214 (15,170)	355,257 (17,744)	414,358 (20,318)	418,293 (20,479)
Unique <i>hkl</i>	83,010 (4,489)	95,751 (4,715)	108,642 (5,344)	100,319 (4,966)
Redundancy	3.4 (3.4)	3.7 (3.8)	3.8 (3.8)	4.2 (4.1)
Completeness (%)	99.8 (99.5)	99.1 (98.4)	99.4 (98.8)	100.0 (100.0)
R <sub>meas</sub>	0.104 (0.678)	0.261 (1.042)	0.078 (0.863)	0.107 (0.678)
CC <sub>1/2</sub>	0.995 (0.750)	0.947 (0.537)	0.998 (0.714)	0.995 (0.769)
I/( $\sigma$ )	8.3 (2.2)	26.6 (4.4)	10.7 (2.0)	7.7 (1.9)
Wilson B (Å <sup>2</sup> )	27.3	27.9	19.1	18.2
R <sub>work</sub> (# <i>hkl</i> )	0.167 (78,804)	0.191 (90,638)	0.189 (103,012)	0.186 (95,349)
R <sub>free</sub> (# <i>hkl</i> )	0.200 (4,183)	0.229 (4,831)	0.230 (5,380)	0.229 (4,951)
B-factors (Å <sup>2</sup> )	32.7	29.9	28.5	24.3
(# atoms)	(7,432)	(7,555)	(7,721)	(7,694)
Ramachandran				
allowed (%)	99.9	100	99.3	99.5
generous (%)	0.1	0	0.7	0.5
disallowed (%)	0	0	0	0
RMSD				
bonds (Å)	0.019	0.020	0.020	0.019
angles (°)	1.89	1.88	1.93	1.86
PDB code	6MXR	6MXS	6MY4	6MY5

the two Fabs (Figure 4(a)), engaging mainly side-chain atoms of charged (Arg), polar (Ser, Asn) and aromatic (Tyr, Trp) residue types, along with a few H-bonds implicating main-chain atoms. This corresponds to  $\sim 200$  Å<sup>2</sup> of buried surface area per direct intermolecular H-bond, which also clearly assigns this interface to the solution homodimer class and not to crystal packing.<sup>25</sup> Additionally, as exchange kinetics are sensitive to changes in the H-bonding network, the two main-chain H-bonds involving backbone amides should be detectable by HX-MS. Protection from exchange was observed in 2 peptides covering T<sup>L</sup>94(NH) (Y<sup>L</sup>87-V<sup>L</sup>114 and Y<sup>L</sup>87-F<sup>L</sup>115), which is involved in H-bonding with S<sup>L</sup>30b(C = O). The H-bond between F<sup>H</sup>98(NH) and Y<sup>H</sup>52(C = O) was not detected due to a gap in sequence coverage. A cavity buried at the dimerization interface contains 13 water molecules and 2 atoms interpreted as chloride ions. These solvent molecules establish further water-mediated H-bond bridges between the Fab monomers (Figure S6). This “wet” region of the interface is in sharp contrast with the rest of the interface that is “dry”, consisting of densely packed solute atoms with absolutely no interfacial water molecules. A total of 24 aromatic residues are found in this “dry” region of the dimer interface (7 Tyr, 3 Trp and 2 Phe from each monomer), engaging in a variety of aromatic ring-stacking modes, hydrophobic contacts and H-bonds (Figure 4(a)). Aliphatic groups are rare at the interface, and are limited to 6 hydrophobic side-chains (2 Ile and 1 Met from each Fab) and 6 methyl groups (from 3 Thr of each Fab). This dimeric Fab interface appears to fit the “core-support-rim” protein binding model of Levy et al.<sup>33</sup> since: 1) the number of amino-acid residues involved in each of the core, support and rim categories are comparable; and 2) the buried surface area of the interface core is significantly larger than those of interface support and rim regions.

### Local conformational changes upon mutation promote aggregation

In addition to the dimerizing conformation, the CDR of the bH1-Fab variant W<sup>H</sup>33-F<sup>H</sup>98-M<sup>H</sup>99 can adopt at least one other low-energy conformation, i.e., the antigen-binding conformation. An overlay between these two conformations is shown in Figure 5. Here, a model is shown of the antigen-bound conformation that was previously generated during the ADAPT affinity maturation campaign using the antigen-bound crystal structure of the parental Fab as a template.<sup>43–45</sup> Relative to the antigen-bound conformation, the most significant changes seen upon dimerization are located in the CDR-H3 loop, which was targeted during affinity maturation and also establishes extensive direct contacts between the Fab monomers. The CDR-H3 conformational change is accompanied by conformational changes in the CDR loops L1 and L3, which have a less intimate involvement in the Fab dimerization interface. The flexibility of the solvent-exposed loop CDR-L1 has been previously reported.<sup>44,45</sup>

The two sets of CDR conformations are specific for their cognate partners. When superimposed onto the antigen-bound structure, the dimerization conformation would clash with VEGF. Similarly, two Fab molecules in the antigen-binding conformation would clash if they approach each other in the observed dimerization mode (Figure S7). A rotation in the CDR-H3 main-chain torsion angles around the Y<sup>H</sup>100a-A<sup>H</sup>100b peptide bond leads to different exposures of side chains F<sup>H</sup>98, M<sup>H</sup>99, F<sup>H</sup>100 and Y<sup>H</sup>100a (Figure 5), making them fitted either for antigen binding or for self-association. The duality of CDR-H3 loop conformations may rely on the fact that both M<sup>H</sup>99 or F<sup>H</sup>100 are capable of strong non-polar packing into an intramolecular hydrophobic pocket, stabilizing either the dimerization or the antigen-binding conformations, respectively. The CDR-



**Figure 5.** Local CDR conformational changes required for self-association. The variable domain of the bH1-Fab variant  $W^{H33}-F^{H98}-M^{H99}$  is shown on the left as an overlay between its dimerization conformation (crystal structure) and antigen-bound conformation (modeled on the antigen-bound structure of the parental bH1-Fab; PDB code 3BDY). The largest conformational changes are seen for three CDR loops as detailed on the right.

H3 antigen-binding conformation might be slightly more stable than the dimerizing conformation, since this same conformation has been previously observed in complex with two different antigens of completely unrelated structures.<sup>44,45</sup>

### Structure-native-aggregation relationship along the affinity maturation path

The ADAPT platform affinity matures an existing antibody by focusing on only a few CDR amino-acids and requires stepwise generation of only a small set of mutants, typically 30–50.<sup>43,50</sup> That self-association can be observed within such a limited set of variants is noteworthy and reinforces recent views on the vulnerability of proteins to self-assembly.<sup>14,51–53</sup> Analysis of the points along this affinity maturation path at which self-association occurs could be instructive and advise as to how readily aggregation is induced by a single mutation. To this end, we reconstituted the ADAPT trail of heavy-chain mutants from the parental bH1-Fab leading to the affinity matured species including the self-associating triple mutants. We also included another affinity-maturing mutation at a fourth position,  $S^{L30bR}$  in the light chain. Table 2 presents experimental data for self-association propensities in this SnAggR set consisting of the parental bH1-Fab and 32 affinity-matured variants with systematic substitutions at 4 CDR positions.

The greatest improvement in the antigen binding affinity resulted from the presence of a hydrophobic amino acid (Met, Trp or Phe) at position  $^H98$ . Secondary selections along the maturation path involved mutations at positions  $^H33$  and  $^H99$ . None of the designed single or double-point mutants of the parental Fab at these three positions show signs of self-association. However, given any combination of designed mutations at two of these sites, designed mutations at the third site triggers self-association. Given hydrophobic residues at  $^H98$  and  $^H99$ , position  $^H33$  becomes a switch that may trigger self-association. Even the conservative substitution  $Y^{H33W}$  in this context is sufficient to initiate self-association. The double mutant with  $Y^{H33W}$  in addition to a hydrophobic residue at  $^H98$  is not prone to self-association, but in this context, a single mutation at position  $^H99$  to a Met residue renders the molecule prone to self-association. The double mutant containing both  $Y^{H33W}$  and  $G^{H99M}$  is not prone to self-association. In this case, it is mutation of position  $^H98$  to a hydrophobic residue that becomes the trigger for self-association. The double-point mutants of the parental Fab can be regarded as single-point mutants with respect to the aggregating triple mutants, and as such underscore the fact that single mutations can become the tipping points toward self-association. The crystal structures of dimeric triple mutants shed light on the self-associating roles of residues at these positions.

As mentioned earlier, the dimerization interface determined here obeys the core-support-rim protein interaction

Table 2. Structure-native-aggregation-relationship (SnAggR) data set of affinity-matured Fabs.

Affinity maturation cycle	Mutated CDR positions <sup>b</sup>										T <sub>m</sub> (°C)	DSC	DLS (10 <sup>-7</sup> cm <sup>2</sup> /s) <sup>d</sup>	dimer % <sup>e</sup>		K <sub>D</sub> dimer (μM)		Res. (Å)	
	H33	H98	H99	L30b	DSF <sup>c</sup>	AUC	AUC	AUC	AUC	ITC				NC					
	Reference																		
Parent	Y	D	G	S	81.1	77.4	7.9	2.8	16000	NB	NC								
Singles	W	D	G	S	80.8		7.8	1.2	1469000										
	Y	M	G	S	80.9		7.8	1.4	179000										
	Y	F	G	S	80.6		7.6	1.0	1290										
	Y	W	G	S		76.7	6.9	2.8	5730										
	Y	D	M	S	79.9		8.1	2.4	220000										
	Y	D	G	R		76.7	7.7	2.6	8940										
Doubles	W	M	G	S	80.6		7.8	1.8	5120000										
	W	F	G	S		76.7	7.4	3.3	3210										
	W	W	G	S		76.3	7.7	1.4	2680										
	Y	M	M	S	80.6		7.8	2.1	19900										
	Y	F	M	S	80.2		7.5	8.1	499										
	Y	W	M	S	79.8		7.5	1.1	725										
	W	D	M	S	80.4		7.8	2.5	35800										
	Y	F	G	R	80.0		7.7	0.6	12800										
Triples	W	M	M	S	79.7	76.4	5.9	86.1	27	8.4	2.04								
	W	F	M	S	79.7	76.2	5.0	93.0	15	4.2	1.95								
	W	W	M	S	79.5	75.9	6.5	78.6	74	45									
	W	M	D	S	80.5	76.7	7.8	1.8	2910										
	W	F	D	S	79.8	76.5	7.8	0.9	15400										
	W	W	D	S	79.6	76.1	7.6	2.7	3240										
	Y	M	D	S	80.1	76.3	7.3	0.3	1210000										
	Y	M	D	R	80.7	77.0	7.8	0.6	185000										
	Y	W	D	R	79.3	76.0	7.6	4.1	4530										
	Y	W	G	R		75.7	6.7	46.7	132	174									
	W	M	G	R	80.4	75.7	7.6	1.6	5610										
	W	F	G	R	80.4	75.7	7.9	1.0	72500										
Quads	W	M	M	R	79.6	75.6	6.1	89.2	13	3.1	1.69								
	W	F	M	R	79.2	75.2	5.8	93.5	13	2.1	1.73								
	W	W	M	R	80.0	75.1	6.1	86.7	36	22									
	W	M	D	R	79.9	76.0	7.7	0.9	3980										
	W	F	D	R	80.0	75.7	8.1	1.1	79100										
	W	W	D	R	79.0	75.2	7.6	1.2	1330										

<sup>a</sup> Antigen binding affinity data taken from ref.<sup>43</sup> except for the Triples and Quads with arginine at position 30b determined in this study.  $\Delta\Delta G_{\text{bind}} = \Delta G_{\text{bind}}(\text{Mutant}) - \Delta G_{\text{bind}}(\text{Parent})$ .

<sup>b</sup> Chothia numbering scheme. Amino-acids mutated relative to the Parent are highlighted.

<sup>c</sup> See Table S4 for standard deviations.

<sup>d</sup> After 14 days of storage at 40°C. Concentration dependence of diffusion coefficients is shown in Figure S11 that also includes data at time-point zero and at 4°C. Diffusion coefficients at tested concentration, as well as calculated self-diffusion coefficients and diffusion interaction parameters are listed in Tables S2.

<sup>e</sup> Determined at 1 mg/mL concentration.


NB: no binding. NC: no crystals.

model.<sup>33</sup> The component interface categories according to this model are shown in Figure 6. An interesting question is whether, starting with this *a priori* determined assembly mode, the interface core propensity (ICP) score associated with the core-support-rim model<sup>33</sup> can explain the self-association trends in the SnAggR dataset. We found that the ICP score calculated only for the interface core residues is not able to correctly classify this data set in terms of self-association propensity (Figure 6, right upper plot). Notably, including interface support residues in the interface propensity calculations improves the classification (Figure 6, right bottom plot). However, some variants remain over-predicted to self-associate even with the combined interface core+support propensity (ICSP) score.

The main culprit for this misclassification is the interface support residue W<sup>H33</sup>. It packs tightly against three aromatic rings from the partner Fab molecule and two intramolecular aromatic residues, each occurring twice per Fab dimer due to the symmetry of the interaction mode (Figure 4(b)). Molecular modeling indicated that a Tyr side-chain would fit less well energetically in this stacking mode, with most of the calculated loss of binding affinity upon W<sup>H33</sup>Y mutation arising from van der Waals contributions (Figure S8). Calculations also indicated that Phe at position <sup>H33</sup> would lead to a further reduction in the tendency for self-association due to the additional loss of an inter-chain H-bond contact

(Figure S8). However, even after including the support category in the propensity calculation, the ICSP score fails to discriminate Y<sup>H33</sup> variants, which are predominantly monomeric, from self-associating W<sup>H33</sup> containing variants (Figure 6). Understandably, the stickiness scale,<sup>33</sup> used for ICSP calculations, is unable to capture the intricate details of short-range interactions established in the complex.

Residue M<sup>H99</sup> also belongs to the support category, having a small surface area buried upon self-association. The crystal structures of the dimeric Fab variants show that the M<sup>H99</sup> side chain does not interact with the partner Fab, but rather with an intramolecular hydrophobic pocket (Figure 4(b)). Introducing the charged Asp side chain or the parental Gly residue at this position completely eliminates self-association, despite the fact that these substitutions appear not to introduce direct steric or electrostatic conflicts in the observed dimerization mode. Inclusion of this support residue in the interface propensity score calculation led to improved correlation with the observed trends in the SnAggR dataset, even though over-prediction of self-association for some of the variants (e.g., bearing G<sup>H99</sup>) still remains. One possible explanation is the conformational freedom of the CDR-H3 loop described earlier, with destabilization of the dimerization conformation by mutating M<sup>H99</sup> to flexible (G<sup>H99</sup>) or charged (D<sup>H99</sup>) residues.



**Figure 6.** Correlations of experimental self-association propensities in the SnAggR dataset with calculated interface propensities based on *a priori* knowledge of the self-assembly mode. Experimental AUC data on dimer fraction is taken from Table 2. Surface categories and interface core propensity (ICP) and interface core+support propensity (ICSP) scores are based on the core-support-rim protein interaction model of Levy et al.<sup>33</sup> Key positions in the SnAggR dataset are delineated by white dotted contours on the surface map on the left.

At position <sup>H</sup>98, mutating Asp to hydrophobic residues converts a completely monomeric Fab (bearing W<sup>H</sup>33 and M<sup>H</sup>99) into a self-associating protein. There appears to be some tolerance at this position for several hydrophobic side chains, as either Met, Phe or Trp lead to Fab dimers, albeit with different self-association strengths that vary in the order Trp < Met < Phe. The crystal structures show that position <sup>H</sup>98 interacts intermolecularly with itself (Figure 4(b)). In the parent molecule, two charged Asp side chains at this position would face one another and may lead to an electrostatic incompatibility in addition to the loss of non-polar packing, thus explaining the absence of self-association of the D<sup>H</sup>98-containing variants. A hydrophobic residue at position <sup>H</sup>98 clearly belongs to the interface core, with interface propensity scores accounting for the SnAggR dataset trends at this position.

The S<sup>L</sup>30bR mutation in the CDR-L1 loop of the light chain appears to further strengthen the affinity of the dimeric Fab variants (Table 2). While the increased presence of hydrophobic residues exposed at the protein surface is typically considered indicative of stronger self-association,<sup>11,17</sup> this mutation demonstrates that other factors may also come into play. In order to understand the structural basis for this behavior, we solved the crystal structures of two Fab variants that incorporate R<sup>L</sup>30b along with W<sup>H</sup>33, M/F<sup>H</sup>98 and M<sup>H</sup>99. These structures show that the introduced Arg side chain of each Fab in the dimer establishes salt-bridge contacts with two Asp side chains, one each from the light and heavy chains of the partner Fab (Figure 4(c)). Interestingly, the <sup>L</sup>30b position changes classification from interface core for the parental Ser to interface rim for the mutant Arg. In the context of R<sup>L</sup>30b, the variant with Arg at position <sup>H</sup>99 showed a reasonable level of dimerization (Table 2), possibly owing to the aliphatic portion of its side chain and to a possible additional intramolecular salt-bridge contact with D<sup>H</sup>101. Notably, this variant introducing two Arg side chains, R<sup>L</sup>30b and R<sup>H</sup>99, partially dimerized in the presence of Y<sup>H</sup>33 that is suboptimal for dimer formation, as discussed earlier.

### Performance of *in silico* methods for aggregation prediction

In contrast with the above analysis based on interface propensity scores, aggregation predictions in typical real-life prospective research have to be made without *a priori* knowledge of the self-association structure. Using the SnAggR data set, we evaluated several popular aggregation prediction methods. The most relevant for our study of self-association between folded monomers are the 3-dimensional (3D)-structure-based methods. In this category, we evaluated 4 approaches: 1) the Spatial Aggregation Propensity (SAP) method based on molecular dynamics (MD) ensembles;<sup>11,16</sup> 2) the SAP-related Developability Index (DI),<sup>17</sup> which was fitted to the present data; 3) the AGGRESCAN3D method in both static and dynamics modes;<sup>54</sup> and 4) the structurally corrected version of the CamSol algorithm.<sup>55</sup> Although sequence-based methods are more suitable for predicting aggregation via the unfolding pathway, including the propensity to form  $\beta$ -amyloid fibrils, we nevertheless tested several of the

many methods available in this category, including AGGRESCAN,<sup>56</sup> CamSol,<sup>55</sup> Waltz,<sup>57</sup> TANGO,<sup>58</sup> and PASTA.<sup>59</sup> We note that this is not an exhaustive list, and interested practitioners may want to use our SnAggR data set for their own evaluation of other preferred methods. Aggregation prediction methods were assessed for their performance in terms of classification of aggregation propensities measured as the fraction of dimer at 1 mg/mL, as well as for more quantitative ranking of aggregation propensities expressed as the logarithm of measured dimer dissociation constant. Performance metrics are given in Table 3, with prediction data for all variants listed in Table S5.

As can be seen, predictions obtained with the tested structure-based methods showed only weak quantitative correlations with experimental aggregation propensities in our systematic SnAggR data set (Table 3, Figure 7(a), Figure S9). The performance for some of these methods increased to a moderate level when employed as classification tools, as indicated by the area under receiver-operating-characteristic (ROC) curve values reaching ~0.8 (ranging from 0.5 for random to 1 for perfect classifications). Notably, these performance levels for both quantitative and qualitative predictions were not superior to those that could be obtained by some sequence-based methods (Table 3, Figure 7(b), Figure S10). Clearly, there was no advantage of including the 3D-based properties that were calculated on unbound structures, i.e., without evaluating possibilities to form self-assembled states. This suggests flexible protein-protein self-docking as a possible way towards improving native-state aggregation predictions.

### Discussion

The notion that hydrophobic surface patches drive protein aggregation is widely accepted and has motivated efforts to predict hot spots for antibody aggregation.<sup>16</sup> In a test case of the predictive method, mutations in the CDR loops were found to meet the selective criteria of reducing aggregation, but sometimes at the expense of antigen interaction.<sup>16</sup> In contrast, our study arises from an affinity maturation campaign in which antibody-antigen interactions were improved

**Table 3.** Performance of various aggregation prediction methods on the SnAggR set.<sup>a</sup>


Prediction Methods <sup>b</sup>	R <sup>2</sup> <sup>c</sup>	$\rho$ <sup>c</sup>	Area under ROC curve <sup>d</sup>
<b>Structure-based</b>			
SAP	0.23	0.48	0.77
DI (Fitted)	0.25	0.50	0.86
CamSol Structure Corrected	0.25	0.49	0.81
AGGRESCAN3D Static	0.00	0.03	0.56
AGGRESCAN3D Dynamic	0.03	0.06	0.55
<b>Sequence-based</b>			
CamSol	0.29	0.54	0.82
AGGRESCAN	0.21	0.49	0.80
TANGO	0.35	0.44	0.82
Waltz	0.00	0.01	0.46
PASTA	0.05	0.25	0.68

<sup>a</sup> Raw data given in Table S5.

<sup>b</sup> See Methods section for details.

<sup>c</sup> Squared Pearson correlation coefficients and Spearman rank-order correlation coefficients are based on dimer LogK<sub>d</sub>( $\mu$ M) data determined by analytical ultracentrifugation and listed in Table 2 and Table S5.

<sup>d</sup> Area under ROC curve (see also Figure 7(b)) values are based on dimer % data determined by analytical ultracentrifugation and listed in Table 2 and Table S5.



**Figure 7.** Performance of various popular *in silico* methods for aggregation predictions on the SnAggr data set without *a priori* knowledge of the self-assembly mode. (a) Scatter plots of experimentally measured self-association propensities expressed as the logarithm of dimer dissociation constant ( $K_d$ ) determined by analytical ultracentrifugation (Table 2) versus aggregation predictions from 3D-structure based methods (see methods section for details). Shown data is generated using the antigen-bound conformation of the CDR (PDB code 3BDY). Alternate predictions generated using the homo-dimeric conformation of the CDR (determined in this work) are provided in Figure S9. (b) Performance of structure and sequence-based methods in terms of classification of aggregation propensities. Receiver operating characteristic (ROC) curves are shown based on the 7 positives (>40% dimer) and 26 negatives (<10% dimer) in the SnAggr data set, with dimer fractions at 1 mg/mL concentration determined experimentally by analytical ultracentrifugation (Table 2). See Table 3 for associated area under ROC curve values. Predicted values for all variants are listed in Table S5.

by mutation of residues in the CDR loops of the antibody.<sup>43</sup> We have shown that selection based on affinity introduced some hydrophobic residues that inadvertently generated

a propensity for symmetric Fab-Fab interactions that lead to aggregation of the antibodies.

Multiple modes of antibody self-association including Fab-Fab, Fab-Fc and Fc-Fc interactions have been documented.<sup>60</sup>

Since affinity maturation mutations studied here are in the CDR, Fc-Fc associations are not expected to play a role in self-association in this case. The strong correlation between mAb aggregation and Fab dimerization (Figures 2 and 3) indicates that the Fab-Fc interactions are also not a major factor.

Our structural data revealed that the antigen-recognition loops (CDR) are directly implicated in the Fab-Fab symmetric interaction interfaces examined here. Large surface areas become buried and the majority of interactions are hydrophobic, although some polar residues and interactions as well as water-mediated contacts are also involved. These “sticky” surface patches formed mainly by aromatic residues do not interact non-specifically, but instead require precise aromatic stacking and are intertwined with short-range polar contacts such as H-bonds and salt-bridges across the interface. A noteworthy finding of this study, however, is that critical residues for self-association do not have to make direct contact across the interface, but instead seem to act as stabilizers of the dimerizing conformation. The best example is M<sup>H99</sup>, which appears to stabilize the CDR-H3 loop in the dimerizing conformation. With a solvent-exposed hydrophobic side-chain located near, but not engaged in, the homodimeric interface, M<sup>H99</sup> is also classified as a support residue according to the core-support-rim architecture of protein-protein interfaces.<sup>30,33</sup> Thus, another class of critical residues for self-association should also be taken into consideration, which consist of interface support residues located in flexible regions (such as antibody CDR-H3 loops) and which are capable of stabilizing self-associating conformations.

The favorable contribution of surface flexibility towards allowing productive conformations for self-association can also benefit from symmetry amplification. Another example revealed here is located in the flexible CDR-L1 loop allowing formation of intermolecular salt-bridges upon introduction of a charged Arg residue at position L<sup>30b</sup>, which is formed twice in the homodimeric Fab. The hinge region connecting the Fab and Fc portions of full-size, homo-dimeric, IgG antibodies is another flexible region that not only enhances binding avidity, but can also modulate the size and shape of the higher-order homooligomers that can be created on the basis of symmetric contacts of each Fab arm, as revealed by electron microscopy images.

Given that low-energy self-assemblies are more likely symmetrical,<sup>37</sup> it is reasonable to expect that self-association by symmetric contacts is not limited to the variants studied here, but is a widespread mechanism of natively-folded aggregation. General support for this hypothesis is found in a recent study showing that symmetry amplifies supramolecular self-assemblies.<sup>14</sup> The “all-or-nothing” behavior observed upon mutation of key positions can be explained on the basis of symmetry amplification of residue contributions to binding affinity, which also contributes to the high occurrence of symmetric homo-oligomeric proteins in nature.<sup>14</sup> Mutating a residue that makes contacts across a symmetric interface, including self-contact, has a doubly disruptive effect relative to mutating a non-symmetric contact. On the basis of symmetry amplification, it is thus not surprising that, in the previously confirmed cases of symmetric self-assemblies, disaggregation was achieved upon mutation of only a few residues.<sup>38–40,42</sup>

Symmetry effects may also contribute to the buffering effect of Asp and other charged residues in protecting against aggregate formation.<sup>14,61</sup> Introduction of Asp at the self-interacting position H<sup>98</sup> of the Fab-Fab interface observed here provides a direct illustration of this buffering phenomenon at the structural level.

The interfaces generated in the affinity maturation exercise have most of the features of known interfaces – buried surface area is significant and comparable to so-called specific interfaces, a cluster of hydrophobic residues contributes significantly to the buried surface, a conformational change contributes to the interaction, complementary surfaces are generated, and few amino acid changes are required to tip the balance toward or away from aggregation. Despite all of those previously observed factors in protein-protein interaction, the sequence and structure-based predictive methods failed to reliably identify the self-aggregating variants. Our structural data is further enriched by aggregation measurements for a set of 33 variants with systematic mutations at key positions of the self-assembly interface. A SnAggR data set with this level of detail, which significantly augments the previous mutagenesis data on structurally characterized Fab self-assembled interfaces,<sup>38–40,42</sup> should be extremely useful for improving computational tools for native aggregation predictions.

Our study demonstrates that improvements must be made with respect to several aspects. Firstly, aggregation predictions with popular structure-based methods, which profile the folded monomer without evaluating possible self-assembled states, did not show quantitative correlations with the data in the present SnAggR set (Figure 7(a) and Figure S9). In fact, there was no clear advantage of including the 3D-structural information, as some of the sequence-based methods showed similar performances (Table 3, Figure 7(b) and Figure S10). This indicates that profiling the monomer structure may not be sufficient, and suggests that further improvements may necessarily require incorporation of structural predictions of protein-protein complexes, i.e., protein docking. Secondly, in the present case, 100-ns MD simulations were unable to achieve conversion of the antigen-bound conformation to the dimerization conformation of the CDR (data not shown), indicating that conformational sampling of the monomeric state *via* room-temperature MD may be insufficient for accurate prediction of self-interactions. Thirdly, short-range polar interactions do contribute to the self-assembly interface and subtle changes, such as residues with very similar hydrophobicities (e.g., Trp, Tyr, Phe), can lead to vastly different aggregation behaviors due to differences in packing and even H-bonding requirements. Exclusive focus on exposed hydrophobicity and coarse electrostatics is thus not sufficient for accurate prediction of aggregation. Self-docking simulation of thoroughly explored ensembles of surface conformations emerges as the appropriate approach to address the required level of complexity observed in this study. Significant efforts have been deployed to improve these computational technologies with some progress currently being made.<sup>62–65</sup> Undoubtedly, protein engineering and antibody optimization will benefit greatly from more accurate predictions of native aggregation. Finally, most computational methods deal with the mechanistic aspects of aggregation and, thus, are meant to reflect the thermodynamics of reversible equilibria such as non-

covalent binding of natively folded monomers (as in this work) or reversible partial unfolding equilibria. However, within the complex pathways of aggregation, post-nucleation events and irreversible unfolding in non-native aggregation are mainly under kinetic control.<sup>4</sup> Experimental measurement of kinetic effects uncoupled from thermodynamically controlled stages of aggregation would be useful for devising and calibrating computational tools intended to also capture kinetic contributions to the overall aggregation process.

## Materials and methods

### Protein production

cDNA for the heavy and light chains of Fab and full-size human IgG1/ $\kappa$  antibody variants of bH1 were ordered from commercial vendors (Thermo-Fisher/Life Technologies Inc., Burlington, ON, Canada; GENEART, Regensburg, Germany). These contained signal peptide sequences, and for Fab variants heavy-chain C-terminal His<sub>8</sub> tags. Productions were carried out by co-transfection of CHO-3E7 cells as described previously,<sup>43</sup> at various scales between 200 mL and 1 L. Transfections were performed at a cell density between  $1.8 \times 10^6$  to  $2.0 \times 10^6$  cells/mL with viability greater than 98%. Cells were distributed in 1.0 L to 2.8 L-shaker flasks and transfected with 1  $\mu$ g of total DNA per 1 mL of production (50% of total DNA contained heavy chain and light chain constructs at ratios of 1:1 (w/w)) using PEI MAX<sup>TM</sup> (Polysciences, Inc., Warrington, PA). The final DNA:PEI MAX<sup>TM</sup> ratio was 1:4 (w/w). Cell cultures were incubated for 24 h on an orbital shaking platform at an agitation rate of 110 rpm at 37°C in a humidified 5% CO<sub>2</sub> atmosphere. Twenty-four hours later, the cultures were fed with Tryptone N1 at 1% w/v final and Valproic acid sodium salt at 0.5 mM final concentration and transferred to 32°C for 6 days. Cell density and cell viability were determined by direct counting of cell samples with a Vi-CELL automated cell counting system (Beckman Coulter Life Sciences, Indianapolis, IN) using the trypan blue dye exclusion method.

### Protein purification

Purifications from cell-culture supernatants were performed by immobilized metal-affinity chromatography for Fab variants and protein-A affinity chromatography for the full-size antibodies (FSAs) having a human IgG1 framework. Fab samples were loaded onto a 1 mL HisTrap<sup>TM</sup> Excel column (GE Healthcare Life Sciences, Uppsala, Sweden) equilibrated in HyClone<sup>TM</sup> Dulbecco's phosphate-buffered saline (DPBS). The column was washed with DPBS and Fabs were eluted with 500 mM imidazole in DPBS. Fractions containing the Fabs were pooled and the imidazole buffer was exchanged against DPBS on PD10 columns (GE Healthcare Life Sciences). For FSAs, purifications of cell-culture supernatants were performed by loading the FSAs onto a 5 mL MabSelect<sup>TM</sup> SuRe<sup>TM</sup> column (GE Healthcare Life Sciences) equilibrated in DPBS. The column was washed with DPBS and FSAs were eluted with 100 mM citrate buffer pH 3.6. The  $3 \times 10$  mL fractions containing FSAs were pooled and the

citrate buffer was exchanged against DPBS on CentriPure P100 columns (EMP Biotech, Berlin, Germany). Purified Fabs and FSAs were aseptically filtered through 0.2  $\mu$ m filters. All affinity purified samples were further purified by preparative SEC on Superdex-200 pg columns (GE Healthcare Life Sciences). Selected peak fractions were concentrated by ultrafiltration using Vivaspin<sup>®</sup> 6 centrifugal concentrators with a membrane molecular weight cut off of 10 kDa (GE Healthcare Life Sciences) at 15°C following the manufacturer's instructions. During the process, the protein concentration was monitored on a NanoDrop<sup>™</sup> 2000 spectrophotometer (Thermo Fisher Scientific, Waltham, MA) using absorbance at 280 nm and the calculated specific extinction coefficient of each variant.

### Differential scanning fluorimetry

Differential scanning fluorimetry (DSF) was used to determine the melting temperatures ( $T_m$ ) of the Fab variants. DSF was carried out in a Rotor-Gene 6000 real-time PCR instrument (Corbett Life Science, Mortlake, NSW, Australia). Samples were diluted in DPBS (GE Healthcare Life Sciences) to a final concentration of 0.33 mg/mL. A total volume of 30  $\mu$ L in 0.2 mL thin wall PCR tubes (Axygen, Oneonta, NY) was used. SYPRO<sup>®</sup> Orange (Life Technologies, Burlington, ON, Canada) was diluted 1,000-fold from the 5,000 $\times$  concentrated stock to the working dye solution in DPBS and 15  $\mu$ L was added to 15  $\mu$ L of sample just prior to the experiment. Thermal denaturation was carried out by increasing the temperature from 30°C to 94°C at a rate of 0.06°C/s. Fluorescence intensity, with excitation at 470 nm and emission at 610 nm, was collected at 1°C intervals and was analyzed with the Rotor Gene 6000 series software v1.7 (Corbett Life Science). The  $T_m$  values were determined from the peak of the first derivative transformation of the raw data. Each sample was measured at least in triplicate and the average  $T_m$  values and standard deviations were reported (Table S4).

### Differential scanning calorimetry

DSC was used to determine the thermal transition midpoint ( $T_m$ ) of selected Fabs and of FSAs. DSC experiments were performed using a VP-Capillary DSC system (Malvern Instruments Ltd, Malvern, UK). Samples were diluted in DPBS (GE Healthcare Life Sciences) to a final concentration of 0.4 mg/mL. Thermal denaturation was carried out under 70 psi of nitrogen pressure by increasing the temperature from 20°C to 100°C at a rate of 60°C/h, with feedback mode/gain set at "low", filtering period of 8 s, prescan time of 3 min. The experiment was run three times with the parental Fab as reference in order to report the precision of the method, but only one time for each other sample. All data were analyzed with Origin 7.0 software (OriginLab Corporation, Northampton, MA). Thermograms were corrected by subtraction of corresponding buffer blank scans and normalized to the protein molar concentration. The  $T_m$  were determined using the automated data processing with the rectangular peak finder algorithm for  $T_m$ .

### Dynamic light scattering

DLS was used to determine the diffusion coefficient of the Fab and FSA variants. Diluted samples were prepared in duplicate in 384-well plate by serial dilution in DPBS (GE Healthcare Life Sciences) to obtain 5 dilution data points ranging from 0.5 to 16 mg/mL. The plates were heat sealed (2.7 s at 170°C, 88 psi) using a PlateLoc thermal microplate sealer (Agilent Technologies, Santa Clara, CA). A second 384-well plate was similarly prepared and incubated for 14 days at 40°C. DLS experiments were performed using a Zetasizer APS system (Malvern Instruments Ltd) equipped with an MRC 150 thermoelectric recirculating chiller (Laird Technologies Inc., Earth City, MO) to thermostat the plate holder at 10°C. The plates were centrifuged for 2 min at 500 × *g* in an Allegra X-14R centrifuge equipped with a SX-4750A rotor (Beckman Coulter, Brea, CA) prior to measurements. The measurements were carried out in triplicate at 25°C after 120 s equilibration time, with automatic measurement duration, number of runs and attenuation factor. All data were analyzed with the Zetasizer software v7.11 (Malvern Instruments Ltd, Malvern, UK), using the Protein Analysis algorithm. After data curation based on the software built-in quality criteria for measurement, the diffusion coefficients were plotted as a function of the sample concentration. Median values of diffusion coefficients ( $D_{\text{median}}$ ) from replicate measurements at each tested concentration were used (Table S1 and Table S2). The self-diffusion coefficient ( $D_0$ ) was determined from the intercept at zero concentration and the diffusion interaction parameter ( $k_D$ ) was calculated from the slope divided by the self-diffusion coefficient. The intercept and slope were obtained through a weighted linear fit attributing a weight proportional to the concentration to account for the higher data variability at lower concentrations (Figure S10). A requirement of at least three data points at a given concentration was imposed for a linear fit to be included.

### Analytical ultracentrifugation

Sedimentation velocity experiments for Fab and FSA variants were performed on a Beckman XL-I analytical ultracentrifuge monitoring absorbance at 280 nm at a protein concentration of 1 mg/mL (~20 μM for Fabs and ~6.8 μM for FSAs) in PBS. Centerpieces were 0.3 cm charcoal-filled epon. Fabs were centrifuged at 48000 rpm and mAbs at 40000 rpm using an 8-hole rotor. The *c(s)* distributions were obtained using SEDFIT,<sup>66</sup> and integrated using GUSI.<sup>67</sup> Binding parameters were obtained by fitting of the data to a monomer-dimer equilibrium model in SEDPHAT.<sup>68</sup> For self-associating variants and their counterparts that differed only by a point mutation, data was collected from 2 different protein lots. For all other Fabs, only a single replicate from a single lot of protein were analyzed. FSAs were analyzed before and after concentration to ~15–20 mg/mL.

### Isothermal titration calorimetry

ITC was performed using a Malvern Microcal Auto VP-ITC (Westborough, MA) to determine the equilibrium dissociation constants,  $K_D$ , for select Fab homodimers. All experiments were carried out using a cell temperature of 25°C, 10 μCal/s reference

power and a stir speed of 1000 rot/min. The Fab variants at concentrations between 100 to 300 μM were titrated into matching buffer in 32 intervals of 1 μL. The resulting isotherms were analyzed using the dissociation model within Origin 7.0 software containing the Microcal macros and the v1.00 dissociation model update.

### Surface plasmon resonance

VEGF antigen binding affinities for the Fab variants bearing S<sup>L</sup>30bR mutation (Table 2) were determined with the same protocol as described for the other Fab variants in this set.<sup>43</sup> Briefly, surface plasmon resonance assays were carried out on a BioRad ProteOn XPR36 instrument (BioRad Laboratories Ltd., Mississauga, ON) at 25°C using PBS running buffer containing 0.05% Tween 20 (Teknova, Hollister, CA) with the addition of 3.4 mM ethylenediaminetetraacetic acid. The GLC sensorchips (Bio-Rad) were prepared by injecting a 1:10 dilution of the standard BioRad sNHS/EDC solutions (Bio-Rad) for 140 s at 100 μL/min, followed by 2.5 mg/mL solution of VEGF antigen (recombinant purified full-length isoform 165 of human VEGF-A) in 10 mM NaOAc pH 4.5 injected at a flow rate of 25 μL/min until approximately 150 resonance units (RUs) of antigen were immobilized. Each Fab variant was injected at 3-fold serial dilutions (top nominal concentrations of 30, 60 or 120 nM) and a buffer blank was simultaneously injected over the blank and VEGF surfaces at 50 μL/min for 120 s with a 900 s dissociation phase. Sensorgrams were aligned and double-referenced using the buffer blank injection and the blank interspots, and the resulting sensorgrams were analyzed using ProteOn Manager software v3.1. Data were fit to the 1:1 binding model to determine the  $k_{\text{on}}$  (s<sup>-1</sup>M<sup>-1</sup>) and  $k_{\text{off}}$  (s<sup>-1</sup>). The binding affinity ( $K_D$ ) was determined from the ratio of  $k_{\text{off}}$ / $k_{\text{on}}$ . Each Fab variant was injected in triplicate.

### Crystallization, data collection, structure determination and refinement

All four bH1 Fab variants were concentrated to 8–10 mg/mL in DPBS buffer before crystallization setup. Crystallization trials were carried out using the microbatch-under-oil method.<sup>69</sup> Both bH1 variants W<sup>H</sup>33-M<sup>H</sup>98-M<sup>H</sup>99 and W<sup>H</sup>33-F<sup>H</sup>98-M<sup>H</sup>99 were crystallized in the condition containing 0.1 M tri-sodium citrate pH 5.6, 11% PEG-4000 and 8% isopropanol, whereas bH1 variants W<sup>H</sup>33-M<sup>H</sup>98-M<sup>H</sup>99-R<sup>L</sup>30b and W<sup>H</sup>33-F<sup>H</sup>98-M<sup>H</sup>99-R<sup>L</sup>30b were obtained in the reservoir solution containing 0.04 M imidazole pH 8.0, 18% PEG-6K. Notably, the parental bH1 Fab did not yield any crystals despite numerous screening trials.

For data collection, the crystals were flash-cooled in liquid nitrogen using the corresponding reservoir solutions supplemented with 15% ethylene glycol as cryo-protectant. Data for the four bH1 Fab variants W<sup>H</sup>33-M<sup>H</sup>98-M<sup>H</sup>99, W<sup>H</sup>33-F<sup>H</sup>98-M<sup>H</sup>99, W<sup>H</sup>33-M<sup>H</sup>98-M<sup>H</sup>99-R<sup>L</sup>30b and W<sup>H</sup>33-F<sup>H</sup>98-M<sup>H</sup>99-R<sup>L</sup>30b were collected at a wavelength of 0.9795 Å to a resolution of 2.04, 1.95, 1.69 and 1.73 Å, respectively, at the CMCF1 beamline, Canadian Light Source. Data sets were indexed and

integrated using iMOSFLM,<sup>70</sup> and scaled using AIMLESS,<sup>71</sup> in the CCP4 suite.<sup>72</sup> Crystals of both W<sup>H33</sup>-M<sup>H98</sup>-M<sup>H99</sup> and W<sup>H33</sup>-F<sup>H98</sup>-M<sup>H99</sup> are isomorphous and belong to the space group P21 with unit cell  $a = 65.7$ ,  $b = 96.0$ ,  $c = 110.1$  Å,  $\beta = 103.7^\circ$ . Both W<sup>H33</sup>-M<sup>H98</sup>-M<sup>H99</sup>-R<sup>L30b</sup> and W<sup>H33</sup>-F<sup>H98</sup>-M<sup>H99</sup>-R<sup>L30b</sup> also crystallized in the P21 space group and have almost identical unit cells ( $a = 81.9$ ,  $b = 65.9$ ,  $c = 93.1$  Å,  $\beta = 99.9^\circ$ ). There are two Fabs (two heavy chains and two light chains) in each asymmetric unit for all these crystals. The structure solution for these Fabs was obtained by molecular replacement through MolRep,<sup>73</sup> using the parental bH1 Fab structure (PDB code 3BDY) as the search model. For each structure, multiple cycles of refinement using REFMAC5,<sup>74</sup> followed by model rebuilding with Coot,<sup>75</sup> were carried out to obtain the final model. All four structures have good stereochemistry as analyzed with PROCHECK.<sup>76</sup> The final models have been deposited into PDB with the accession codes 6MXR, 6MXS, 6MY4 and 6MY5. Data collection and refinement statistics are shown in Table 1.

### Hydrogen-exchange mass spectrometry

Stock solutions (5 mg/mL) were diluted to 40  $\mu$ M with 1  $\times$  PBS (pH 7.9). Deuteration was initiated by a two-fold dilution with 1  $\times$  PBS in 90% D<sub>2</sub>O (final D<sub>2</sub>O 45%). Labeling was performed for 0.05, 0.16, 1, 15 and 60 min at 10°C. The reaction was quenched by a 30-fold dilution with 8 M urea/1 M tris(2-carboxyethyl) phosphine hydrochloride (TCEP) in 100 mM Gly-HCl (pH 3.0), and allowed to reduce for 3 min at 10°C. Samples (0.66  $\mu$ M, 13 pmol) were manually injected into a 19  $\mu$ L loop (13 pmol). Fabs were digested with a Poroszyme immobilized pepsin cartridge (30  $\times$  2.1 mm) at 0°C for 3.5 min (100  $\mu$ L/min, 3% acetonitrile (ACN), 0.1% formic acid (FA), pH 2.7).<sup>77</sup> The resulting peptides were trapped on a C18 PepMap100  $\mu$ precolumn (5  $\mu$ m, 100 Å, 0.3  $\times$  5 mm) and eluted with an 8-min gradient (15–35%) of mobile phase B (97% ACN, 0.1% FA, pH 2.7) at 5  $\mu$ L/min and separated with a BioBasic-18 analytical column (5  $\mu$ m, 300 Å, 50  $\times$  0.32 mm). All columns were purchased from Thermo Fisher Scientific (Waltham, MA). The pepsin column was washed following a procedure proposed previously.<sup>78</sup> Data was collected with a Waters Ultima Q-TOF scanning from 300–1500  $m/z$  coupled to Agilent 1100/1260 LC pumps. Coverage maps were generated with iterative data dependent MS/MS acquisitions with both fully reduced (8 M Urea/1 M TCEP at room temperature overnight) and “HX-compatible” samples.<sup>79</sup> Fitted deuteration was calculated using MS Studio. Data was collected in triplicate, and only peptides found in all replicates were included in the analysis. A 95% confidence interval of  $\pm 0.2$  and 0.5 Da were calculated for individual and summed time points, respectively. Statistical significance was assigned based previously proposed criteria.<sup>80</sup>

### Transmission electron microscopy

Carbon-coated grids (400 mesh, Ted Pella, Redding, CA) were glow discharged at 25 mA for 40 s. Antibody samples were diluted in DPBS to 10  $\mu$ g/ml and 3  $\mu$ L was spotted onto the glow discharged grid and allowed to absorb for 1 min. Excess solution was blotted and the grids were then washed with

MilliQ-H<sub>2</sub>O with blotting on Whatman paper between wash steps, stained with 2% (w/v) uranyl acetate, blotted, air dried and imaged on a Hitachi 7500 TEM (Itachi, Ltd., Tokyo, Japan) operating at 80 keV. Particle selection, alignment and sorting for reference-free 2D class averaging was done in RELION 2.1.0 (MRC Laboratory for Molecular Biology, Cambridge, UK).

### In silico predictions of aggregation propensities

SAP and DI were calculated as previously described<sup>16,17</sup> on the entire panel of variants. SAP calculations were performed using a sphere of radius 5 Å. Surfaces were generated using an in-house molecular surface generator. The different functional forms of DI were evaluated, and only the form with best predictive power is reported. DI was calibrated against sedimentation velocity AUC measurements (% dimer) using least-square linear regression in R.<sup>81</sup> Median SAP values were based on conformational ensembles generated by MD, whereas DI values were calculated on static structures. The structure of the Fab of the variants was modeled using as template X-ray structures either: 1) the antigen-bound conformation (parental bH1-Fab, PDB code 3BDY), or 2) the dimerizing conformation (bH1-Fab variant W<sup>H33</sup>-F<sup>H98</sup>-M<sup>H99</sup>-R<sup>L30b</sup>, chains H and L, this study). The residues were mutated with the PyMOL program (Schrödinger, Inc., New York, NY) and the first rotameric state with steric clearance was selected. The backbone atoms were kept in their initial state. Conformational ensembles of the variants were then generated *via* MD simulations. Initial structures for MD were prepared first with the Sybyl program (Tripos, Inc., St. Louis, MO) to add missing atoms and then using the tleap program from AMBER 16 software.<sup>82</sup> The AMBER FF99SB force field,<sup>83</sup> was used to perform the MD simulations. Each system was solvated in a truncated octahedron TIP3P water box.<sup>84</sup> The distance between the wall of the box and the closest atom of the solute was 12.0 Å. Counterions (Na<sup>+</sup>, Cl<sup>-</sup>) were added up to a final concentration of 0.1 M to maintain electroneutrality of the systems. Each system was minimized first, applying harmonic restraints with force constants of 1 kcal/(mol Å<sup>2</sup>) to all solute atoms, followed by heating from 10 to 150 K for 30 ps in the canonical ensemble (NVT) and from 150 to 300 K for 100 ps in the isothermal–isobaric ensemble (NPT). Each system was equilibrated to adjust the solvent density under 1 atm pressure through 1 ns of NPT simulation. A 30-ns NVT production run was obtained with snapshots collected every 10 ps. For all simulations, a 2-fs time step and an 8-Å nonbonded cutoff were used.

The sequence-based prediction method AGGRESCAN<sup>56</sup> was run through the web server (<http://bioinf.uab.es/aggreSCAN/>) and the  $a^4$ vSS score (window-averaged aggregation propensities values summed over the entire sequence) normalized for 100 residues (Na<sup>4</sup>vSS) was reported. The structure-based AGGRESCAN3D<sup>54</sup> method was run via the AggreSCAN3D 2.0 web server (<http://biocomp.chem.uw.edu.pl/A3D2/>) using a sphere of 10 Å separately for static and dynamic modes.

The CamSol<sup>55</sup> method was run in the sequence-based (intrinsic) mode via the web server (<http://www-mvsoftware.ch.cam.ac.uk/index.php/camsolintrinsic>) and in the structurally corrected mode (<http://www-mvsoftware.ch.cam.ac.uk/index.php/camsol>)

strucorr) using a patch radius of 10 Å. Aggregation propensity scores were obtained by multiplying by a factor of  $-1$  the intrinsic or the structurally corrected scores.

The sequence-based method Waltz<sup>57</sup> was run through its web server (<http://waltz.switchlab.org/>) in the Best Overall Performance mode (92% specificity threshold). The aggregation propensity score was obtained from the sum of scores over all independent aggregation-prone sequence clusters identified.

The sequence-based method TANGO<sup>58</sup> was run locally using the standalone executable version 2.3.1 obtained from the web server (<http://tango.crg.es>) with default settings at 298.15 K. The reported aggregation propensity score was used.

The sequence-based method PASTA<sup>59</sup> was run at the PASTA 2.0 web server (<http://protein.bio.unipd.it/pasta2/>) using the parameters corresponding to 90% specificity threshold. The aggregation propensity score was obtained from the sum of energies for the top 22 self-aggregating pairings, and then multiplied by a factor of  $-1$ .

Contributions to changes in binding affinity upon mutation at position H33 were estimated with the Solvated Interaction Energy (SIE) function.<sup>64,85</sup> The structure of mutants Y<sup>H33</sup> and F<sup>H33</sup> were predicted and modelled in PyMOL from the crystal structure of the dimer Fab with W<sup>H33</sup>. Dihedral angles  $\chi_1$  and  $\chi_2$  were manually fixed to those from the template structure. A constrained minimization protocol was applied such that: 1) the hydrogen atoms were free to move, 2) the side-chains of the mutated residue and residues in close proximity (5 Å) were imposed a 1.0 kcal/(mol Å<sup>2</sup>) harmonic restraint, and 3) all other atoms were imposed a 10.0 kcal/mol kcal/(mol Å<sup>2</sup>) harmonic restraint.

Interface propensity scores (ICP, ICSP) were calculated based on SASA generated with an in-house program and amino-acid propensities summed over the interface core and support surfaces as described previously.<sup>33</sup> The Fab dimer crystal structure of bH1 variant W<sup>H33</sup>-F<sup>H98</sup>-M<sup>H99</sup> was used as reference variant, as well as for the assignment of core, support and rim residues. The scores of the rest of the variants were calculated relative to the score of the reference variant.

## Acknowledgments

We are grateful for the technical assistance of Marie-France Goneau, Ian Schoenhofen, Melanie Duchesne, Mark Sasserville, Andrea Accl, Aline Do-Minh, Nathalie Bousquet-Gagnon, Louise Thauvette, Hongtao Qi, Roger Maurice, Erin Weston, The-Minh Tu, Marie-Josée Massariol, Nancy Remillard, Alina Burlacu and Guillaume Arthus-Cartier. We thank Allan Matte for useful discussions and support. Crystallographic research described in this paper was partially supported by the Natural Sciences and Engineering Research Council (NSERC) of Canada discovery grant 436202 to R.S. and was performed using beamline 08ID-1 at the Canadian Light Source, which is supported by the Canada Foundation for Innovation, NSERC, the University of Saskatchewan, the Government of Saskatchewan, Western Economic Diversification Canada, the National Research Council Canada, and the Canadian Institutes of Health Research.






## Funding

This work was supported by the Natural Sciences and Engineering Research Council of Canada [436202].

## Author contributions

J.D.S. and T.S. designed research; J.D.S., M.E.P., F.G., L.P.G., J.B., M.S.M., J.S. and C.B. performed research; C.D. and H.H. contributed new analytic tools; J.D.S., M.E.P., F.G., L.P.G., J.B., M.S.M., J.S., J.K., E.O.P., R.S. and T.S. analyzed data; J.D.S., M.E.P., F.G., L.P.G., J.B., J.S., R.S. and T.S. wrote the paper and all authors contributed to manuscript editing.

## ORCID

Marie-Ève Picard  <http://orcid.org/0000-0002-4316-2830>  
Francis Gaudreault  <http://orcid.org/0000-0001-5656-3118>  
Louis-Patrick Gagnon  <http://orcid.org/0000-0001-7417-866X>  
Jason Baardsnes  <http://orcid.org/0000-0001-5587-3034>  
Mahder S. Manenda  <http://orcid.org/0000-0002-8640-5556>  
Joey Sheff  <http://orcid.org/0000-0002-9269-6863>  
Christophe Deprez  <http://orcid.org/0000-0003-2687-3144>  
Rong Shi  <http://orcid.org/0000-0001-8494-900X>  
Traian Sulea  <http://orcid.org/0000-0001-5301-8261>

## References

1. Kumar S, Robins RH, Buck PM, Hickling TP, Thangakani AM, Li LS, Singh SK, Gromiha MM. Biopharmaceutical informatics: applications of computation in biologic drug development. In: Kumar S, Singh SK, editors. Developability of biotherapeutics: computational approaches. Boca Raton (FL): CRC Press; 2016. p. 3–34.
2. Agrawal NJ, Chennamsetty N. Understanding, predicting, and mitigating the impact of posttranslational physicochemical modifications, including aggregation, on the stability of biopharmaceutical drug products. In: Kumar S, Singh SK, editors. Developability of biotherapeutics: computational approaches. Boca Raton (FL): CRC Press; 2016. p. 61–84.
3. Chaudhri A. Challenges in high-concentration biopharmaceutical drug delivery: a modeling perspective. In: Kumar S, Singh SK, editors. Developability of biotherapeutics: computational approaches. Boca Raton (FL): CRC Press; 2016. p. 135–52.
4. Roberts CJ. Therapeutic protein aggregation: mechanisms, design, and control. Trends Biotechnol. 2014;32:372–80. doi:10.1016/j.tibtech.2014.05.005.
5. Tomar DS, Kumar S, Singh SK, Goswami S, Li L. Molecular basis of high viscosity in concentrated antibody solutions: strategies for high concentration drug product development. MABS. 2016;8:216–28. doi:10.1080/19420862.2015.1128606.
6. Lee CC, Perchiacca JM, Tessier PM. Toward aggregation-resistant antibodies by design. Trends Biotechnol. 2013;31:612–20. doi:10.1016/j.tibtech.2013.07.002.
7. Lowe D, Dudgeon K, Rouet R, Schofield P, Jeremius L, Christ D. Aggregation, stability, and formulation of human antibody therapeutics. Adv Protein Chem Struct Biol. 2011;84:41–61. doi:10.1016/B978-0-12-386483-3.00004-5.
8. Zurdo J, Arnell A, Obrezanova O, Smith N, Gallagher TRA, de la Cuesta RG, Locwin B. Developability assessment workflows to de-risk biopharmaceutical development. In: Kumar S, Singh SK, editors. Developability of biotherapeutics: computational approaches. Boca Raton (FL): CRC Press; 2016. p. 221–89.
9. Tiller KE, Tessier PM. Advances in antibody design. Annu Rev Biomed Eng. 2015;17:191–216. doi:10.1146/annurev-bioeng-071114-040733.
10. Meric G, Robinson AS, Roberts CJ. Driving forces for nonnative protein aggregation and approaches to predict aggregation-prone regions. Annu Rev Chem Biomol Eng. 2017;8:139–59. doi:10.1146/annurev-chembioeng-060816-101404.
11. Chennamsetty N, Voynov V, Kayser V, Helk B, Trout BL. Prediction of aggregation prone regions of therapeutic proteins. J Phys Chem B. 2010;114:6614–24. doi:10.1021/jp911706q.
12. Agrawal NJ, Kumar S, Wang X, Helk B, Singh SK, Trout BL. Aggregation in protein-based biotherapeutics: computational

- studies and tools to identify aggregation-prone regions. *J Pharm Sci.* 2011;100:5081–95. doi:10.1002/jps.22705.
13. Clarkson BR, Schon A, Freire E. Conformational stability and self-association equilibrium in biologics. *Drug Discov Today.* 2016;21:342–47. doi:10.1016/j.drudis.2015.11.007.
  14. Garcia-Seisdedos H, Empereur-Mot C, Elad N, Levy ED. Proteins evolve on the edge of supramolecular self-assembly. *Nature.* 2017;548:244–47. doi:10.1038/nature23320.
  15. Liu J, Nguyen MD, Andya JD, Shire SJ. Reversible self-association increases the viscosity of a concentrated monoclonal antibody in aqueous solution. *J Pharm Sci.* 2005;94:1928–40. doi:10.1002/jps.20347.
  16. Chennamsetty N, Voynov V, Kayser V, Helk B, Trout BL. Design of therapeutic proteins with enhanced stability. *Proc Natl Acad Sci USA.* 2009;106:11937–42. doi:10.1073/pnas.0904191106.
  17. Lauer TM, Agrawal NJ, Chennamsetty N, Egodage K, Helk B, Trout BL. Developability index: a rapid in silico tool for the screening of antibody aggregation propensity. *J Pharm Sci.* 2012;101:102–15. doi:10.1002/jps.22758.
  18. Pindrus M, Shire SJ, Kelley RF, Demeule B, Wong R, Xu Y, Yadav S. Solubility challenges in high concentration monoclonal antibody formulations: relationship with amino acid sequence and intermolecular interactions. *Mol Pharm.* 2015;12:3896–907. doi:10.1021/acs.molpharmaceut.5b00336.
  19. Chaudhri A, Zarraga IE, Yadav S, Patapoff TW, Shire SJ, Voth GA. The role of amino acid sequence in the self-association of therapeutic monoclonal antibodies: insights from coarse-grained modeling. *J Phys Chem B.* 2013;117:1269–79. doi:10.1021/jp3108396.
  20. Buck PM, Chaudhri A, Kumar S, Singh SK. Highly viscous antibody solutions are a consequence of network formation caused by domain-domain electrostatic complementarities: insights from coarse-grained simulations. *Mol Pharm.* 2015;12:127–39. doi:10.1021/mp500485w.
  21. Sharma VK, Patapoff TW, Kabakoff B, Pai S, Hilario E, Zhang B, Li C, Borisov O, Kelley RF, Chorny I, et al. In silico selection of therapeutic antibodies for development: viscosity, clearance, and chemical stability. *Proc Natl Acad Sci U S A.* 2014;111:18601–06. doi:10.1073/pnas.1421779112.
  22. Brader ML, Baker EN, Dunn MF, Laue TM, Carpenter JF. Using X-ray crystallography to simplify and accelerate biologics drug development. *J Pharm Sci.* 2017;106:477–94. doi:10.1016/j.xphs.2016.10.017.
  23. Venien-Bryan C, Li Z, Vuillard L, Boutin JA. Cryo-electron microscopy and X-ray crystallography: complementary approaches to structural biology and drug discovery. *Acta Crystallogr F Struct Biol Commun.* 2017;73:174–83. doi:10.1107/S2053230X17003740.
  24. Bahadur RP, Chakrabarti P, Rodier F, Janin J. A dissection of specific and non-specific protein-protein interfaces. *J Mol Biol.* 2004;336:943–55.
  25. Bahadur RP, Zacharias M. The interface of protein-protein complexes: analysis of contacts and prediction of interactions. *Cell Mol Life Sci.* 2008;65:1059–72. doi:10.1007/s00018-007-7451-x.
  26. Lo Conte L, Chothia C, Janin J. The atomic structure of protein-protein recognition sites. *J Mol Biol.* 1999;285:2177–98.
  27. Ponstingl H, Henrick K, Thornton JM. Discriminating between homodimeric and monomeric proteins in the crystalline state. *Proteins.* 2000;41:47–57.
  28. Chen J, Sawyer N, Regan L. Protein-protein interactions: general trends in the relationship between binding affinity and interfacial buried surface area. *Protein Sci.* 2013;22:510–15. doi:10.1002/pro.2230.
  29. Luo J, Liu Z, Guo Y, Li M. A structural dissection of large protein-protein crystal packing contacts. *Sci Rep.* 2015;5:14214. doi:10.1038/srep14214.
  30. Chakrabarti P, Janin J. Dissecting protein-protein recognition sites. *Proteins.* 2002;47:334–43.
  31. Janin J. Wet and dry interfaces: the role of solvent in protein-protein and protein-DNA recognition. *Structure.* 1999;7:R277–279.
  32. Rodier F, Bahadur RP, Chakrabarti P, Janin J. Hydration of protein-protein interfaces. *Proteins.* 2005;60:36–45. doi:10.1002/prot.20478.
  33. Levy ED. A simple definition of structural regions in proteins and its use in analyzing interface evolution. *J Mol Biol.* 2010;403:660–70. doi:10.1016/j.jmb.2010.09.028.
  34. Monod J, Wyman J, Changeux JP. On the nature of allosteric transitions: a plausible model. *J Mol Biol.* 1965;12:88–118.
  35. Edelstein SJ. Extensions of the allosteric model for haemoglobin. *Nature.* 1971;230:224–27. doi:10.1038/230224a0.
  36. Grueninger D, Treiber N, Ziegler MO, Koetter JW, Schulze MS, Schulz GE. Designed protein-protein association. *Science.* 2008;319:206–09. doi:10.1126/science.1150421.
  37. Andre I, Strauss CE, Kaplan DB, Bradley P, Baker D. Emergence of symmetry in homooligomeric biological assemblies. *Proc Natl Acad Sci USA.* 2008;105:16148–52. doi:10.1073/pnas.0807576105.
  38. Pepinsky RB, Silvan L, Berkowitz SA, Farrington G, Lugovskoy A, Walus L, Eldredge J, Capili A, Mi S, Graff C, et al. Improving the solubility of anti-LINGO-1 monoclonal antibody Li33 by isotype switching and targeted mutagenesis. *Protein Sci.* 2010;19:954–66. doi:10.1002/pro.372.
  39. Chow CK, Allan BW, Chai Q, Atwell S, Lu J. Therapeutic antibody engineering to improve viscosity and phase separation guided by crystal structure. *Mol Pharm.* 2016;13:915–23. doi:10.1021/acs.molpharmaceut.5b00817.
  40. Bethea D, Wu SJ, Luo J, Hyun L, Lacy ER, Teplyakov A, Jacobs SA, O'Neil KT, Gilliland GL, Feng Y. Mechanisms of self-association of a human monoclonal antibody CNTO607. *Protein Eng Des Sel.* 2012;25:531–37. doi:10.1093/protein/gz047.
  41. Teplyakov A, Obmolova G, Wu SJ, Luo J, Kang J, O'Neil K, Gilliland GL. Epitope mapping of anti-interleukin-13 neutralizing antibody CNTO607. *J Mol Biol.* 2009;389:115–23. doi:10.1016/j.jmb.2009.03.076.
  42. Wu SJ, Luo J, O'Neil KT, Kang J, Lacy ER, Canziani G, Baker A, Huang M, Tang QM, Raju TS, et al. Structure-based engineering of a monoclonal antibody for improved solubility. *Protein Eng Des Sel.* 2010;23:643–51. doi:10.1093/protein/gzq037.
  43. Vivcharuk V, Baardsnes J, Deprez C, Sulea T, Jaramillo M, Corbeil CR, Mullick A, Magoon J, Marcil A, Durocher Y, et al. Assisted design of antibody and protein therapeutics (ADAPT). *PLoS One.* 2017;12:e0181490. doi:10.1371/journal.pone.0181490.
  44. Bostrom J, Yu SF, Kan D, Appleton BA, Lee CV, Billeci K, Man W, Peale F, Ross S, Wiesmann C, et al. Variants of the antibody Herceptin that interact with HER2 and VEGF at the antigen binding site. *Science.* 2009;323:1610–14. doi:10.1126/science.1165480.
  45. Bostrom J, Haber L, Koenig P, Kelley RF, Fuh G. High affinity antigen recognition of the dual specific variants of herceptin is entropy-driven in spite of structural plasticity. *PLoS One.* 2011;6:e17887. doi:10.1371/journal.pone.0017887.
  46. Pokric B, Pucar Z. The two-cross immunodiffusion technique: diffusion coefficients and precipitating titers of IgG in human serum and rabbit serum antibodies. *Anal Biochem.* 1979;93:103–14.
  47. Jossang T, Feder J, Rosenqvist E. Photon correlation spectroscopy of human IgG. *J Protein Chem.* 1988;7:165–71.
  48. Dobson CL, Devine PW, Phillips JJ, Higazi DR, Lloyd C, Popovic B, Arnold J, Buchanan A, Lewis A, Goodman J, et al. Engineering the surface properties of a human monoclonal antibody prevents self-association and rapid clearance in vivo. *Sci Rep.* 2016;6:38644. doi:10.1038/srep38644.
  49. Arora J, Hickey JM, Majumdar R, Esfandiary R, Bishop SM, Samra HS, Middaugh CR, Weis DD, Volkin DB. Hydrogen exchange mass spectrometry reveals protein interfaces and distant dynamic coupling effects during the reversible self-association of an IgG1 monoclonal antibody. *MAbs.* 2015;7:525–39. doi:10.1080/19420862.2015.1029217.
  50. Sulea T, Hussack G, Ryan S, Tanha J, Purisima EO. Application of assisted design of antibody and protein therapeutics (ADAPT) improves efficacy of a *Clostridium difficile* toxin A

- single-domain antibody. *Sci Rep.* 2018;8:2260. doi:10.1038/s41598-018-20599-4.
51. Kumar S, Roffi K, Tomar DS, Cirelli D, Luksha N, Meyer D, Mitchell J, Allen MJ, Li L. Rational optimization of a monoclonal antibody for simultaneous improvements in its solution properties and biological activity. *Protein Eng Des Sel.* 2018;31:313–25. doi:10.1093/protein/gzy020.
  52. Geoghegan JC, Fleming R, Damschroder M, Bishop SM, Sathish HA, Esfandiary R. Mitigation of reversible self-association and viscosity in a human IgG1 monoclonal antibody by rational, structure-guided Fv engineering. *MAbs.* 2016;8:941–50. doi:10.1080/19420862.2016.1171444.
  53. Shan L, Mody N, Sormani P, Rosenthal KL, Damschroder MM, Esfandiary R. Developability assessment of engineered monoclonal antibody variants with a complex self-association behavior using complementary analytical and in silico tools. *Mol Pharm.* 2018;15:5697–710. doi:10.1021/acs.molpharmaceut.8b00867.
  54. Zambrano R, Jamroz M, Szczasiuk A, Pujols J, Kmiecik S, Ventura S. AGGRESCAN3D (A3D): server for prediction of aggregation properties of protein structures. *Nucleic Acids Res.* 2015;43:W306–313. doi:10.1093/nar/gkv359.
  55. Sormani P, Aprile FA, Vendruscolo M. The CamSol method of rational design of protein mutants with enhanced solubility. *J Mol Biol.* 2015;427:478–90. doi:10.1016/j.jmb.2014.09.026.
  56. Conchillo-Sole O, de Groot NS, Aviles FX, Vendrell J, Daura X, Ventura S. AGGRESCAN: a server for the prediction and evaluation of “hot spots” of aggregation in polypeptides. *BMC Bioinformatics.* 2007;8:65. doi:10.1186/1471-2105-8-65.
  57. Oliveberg M. Waltz, an exciting new move in amyloid prediction. *Nat Methods.* 2010;7:187–88. doi:10.1038/nmeth0310-187.
  58. Fernandez-Escamilla AM, Rousseau F, Schymkowitz J, Serrano L. Prediction of sequence-dependent and mutational effects on the aggregation of peptides and proteins. *Nat Biotechnol.* 2004;22:1302–06. doi:10.1038/nbt1012.
  59. Walsh I, Seno F, Tosatto SC, Trovato A. PASTA 2.0: an improved server for protein aggregation prediction. *Nucleic Acids Res.* 2014;42:W301–307. doi:10.1093/nar/gku399.
  60. Plath F, Ringler P, Graff-Meyer A, Stahlberg H, Lauer ME, Rufer AC, Graewert MA, Svergun D, Gellermann G, Finkler C, et al. Characterization of mAb dimers reveals predominant dimer forms common in therapeutic mAbs. *MAbs.* 2016;8:928–40. doi:10.1080/19420862.2016.1168960.
  61. Dudgeon K, Rouet R, Kokmeijer I, Schofield P, Stolp J, Langley D, Stock D, Christ D. General strategy for the generation of human antibody variable domains with increased aggregation resistance. *Proc Natl Acad Sci USA.* 2012;109:10879–84. doi:10.1073/pnas.1202866109.
  62. Huang SY. Exploring the potential of global protein-protein docking: an overview and critical assessment of current programs for automatic ab initio docking. *Drug Discov Today.* 2015;20:969–77. doi:10.1016/j.drudis.2015.03.007.
  63. Weitzner BD, Jeliakov JR, Lyskov S, Marze N, Kuroda D, Frick R, Adolf-Bryfogle J, Biswas N, Dunbrack RL Jr., Gray JJ. Modeling and docking of antibody structures with Rosetta. *Nat Protoc.* 2017;12:401–16. doi:10.1038/nprot.2016.180.
  64. Sulea T, Vivcharuk V, Corbeil CR, Deprez C, Purisima EO. Assessment of solvated interaction energy function for ranking antibody-antigen binding affinities. *J Chem Inf Model.* 2016;56:1292–303. doi:10.1021/acs.jcim.6b00043.
  65. Hogues H, Gaudreault F, Corbeil CR, Deprez C, Sulea T, Purisima EO. ProPOSE: direct exhaustive protein-protein docking with side chain flexibility. *J Chem Theory Comput.* 2018;14:4938–47. doi:10.1021/acs.jctc.8b00225.
  66. Schuck P. Size-distribution analysis of macromolecules by sedimentation velocity ultracentrifugation and lamm equation modeling. *Biophys J.* 2000;78:1606–19. doi:10.1016/S0006-3495(00)76713-0.
  67. Brautigam CA. Calculations and publication-quality illustrations for analytical ultracentrifugation data. *Methods Enzymol.* 2015;562:109–33. doi:10.1016/bs.mie.2015.05.001.
  68. Schuck P. On the analysis of protein self-association by sedimentation velocity analytical ultracentrifugation. *Anal Biochem.* 2003;320:104–24.
  69. Chayen NE, Shaw Stewart PD, Blow DM. Microbatch crystallization under oil - a new technique allowing many small-volume crystallization trials. *J Cryst Growth.* 1992;122:176–80. doi:10.1016/0022-0248(92)90241-A.
  70. Battye TG, Kontogiannis L, Johnson O, Powell HR, Leslie AG. IMOSFLM: a new graphical interface for diffraction-image processing with MOSFLM. *Acta Crystallogr D Biol Crystallogr.* 2011;67:271–81. doi:10.1107/S0907444910048675.
  71. Evans PR, Murshudov GN. How good are my data and what is the resolution? *Acta Crystallogr D Biol Crystallogr.* 2013;69:1204–14. doi:10.1107/S0907444913000061.
  72. Winn MD, Ballard CC, Cowtan KD, Dodson EJ, Emsley P, Evans PR, Keegan RM, Krissinel EB, Leslie AG, McCoy A, et al. Overview of the CCP4 suite and current developments. *Acta Crystallogr D Biol Crystallogr.* 2011;67:235–42. doi:10.1107/S0907444910045749.
  73. Vagin A, Teplyakov A. Molecular replacement with MOLREP. *Acta Crystallogr D Biol Crystallogr.* 2010;66:22–25. doi:10.1107/S0907444909042589.
  74. Murshudov GN, Vagin AA, Dodson EJ. Refinement of macromolecular structures by the maximum-likelihood method. *Acta Crystallogr D Biol Crystallogr.* 1997;53:240–55. doi:10.1107/S0907444996012255.
  75. Emsley P, Cowtan K. Coot: model-building tools for molecular graphics. *Acta Crystallogr D Biol Crystallogr.* 2004;60:2126–32. doi:10.1107/S0907444904019158.
  76. Laskowski RA, MacArthur MW, Moss DS, Thornton JM. PROCHECK: a program to check the stereochemical quality of protein structures. *J Appl Cryst.* 1993;26:283–91. doi:10.1107/S002188982009944.
  77. Wang L, Pan H, Smith DL. Hydrogen exchange-mass spectrometry: optimization of digestion conditions. *Mol Cell Proteomics.* 2002;1:132–38. doi:10.1074/mcp.m100009-mcp200.
  78. Majumdar R, Manikwar P, Hickey JM, Arora J, Middaugh CR, Volkin DB, Weis DD. Minimizing carry-over in an online pepsin digestion system used for the H/D exchange mass spectrometric analysis of an IgG1 monoclonal antibody. *J Am Soc Mass Spectrom.* 2012;23:2140–48. doi:10.1007/s13361-012-0485-9.
  79. Hamuro Y, Coales SJ. Optimization of feasibility stage for hydrogen/deuterium exchange mass spectrometry. *J Am Soc Mass Spectrom.* 2018;29:623–29. doi:10.1007/s13361-017-1860-3.
  80. Houde D, Berkowitz SA, Engen JR. The utility of hydrogen/deuterium exchange mass spectrometry in biopharmaceutical comparability studies. *J Pharm Sci.* 2011;100:2071–86. doi:10.1002/jps.22432.
  81. R Development Core Team. R: a language and environment for statistical computing. Vienna (Austria): The R Foundation for Statistical Computing; 2011.
  82. Case DA, Cheatham TE, Darden T, Gohlke H, Luo R, Merz KM, Onufriev A, Simmerling C, Wang B, Woods RJ. The Amber biomolecular simulation programs. *J Comput Chem.* 2005;26:1668–88. doi:10.1002/jcc.20290.
  83. Hornak V, Abel R, Okur A, Strockbine B, Roitberg A, Simmerling C. Comparison of multiple Amber force fields and development of improved protein backbone parameters. *Proteins.* 2006;65:712–25. doi:10.1002/prot.21123.
  84. Jorgensen WL, Chandrasekhar J, Madura JD, Impey RW, Klein ML. Comparison of simple potential functions for simulating liquid water. *J Chem Phys.* 1983;79:926–35. doi:10.1063/1.445869.
  85. Sulea T, Purisima EO. The solvated interaction energy method for scoring binding affinities. *Methods Mol Biol.* 2012;819:295–303. doi:10.1007/978-1-61779-465-0\_19.