



Supplementary Fig. 3: Internalization of eGFP-DODMA-LNP (eGFP-LNP) is reduced when clathrin-mediated endocytic pathway is inhibited. BMDC were treated with 20 μ M chlorpromazine (CP, Sigma-Aldrich) for 1 hr at 37°C prior to transfection with 3 μ g ApoC3-eGFP-LNP for 24 hr followed by flow cytometry to determine % GFP positive cells. Data is plotted relative to eGFP-LNP. n=3. $P < 0.01$ (*) is relative to eGFP-LNP.