




Sample ID: 18-19159-022
Report ID: R_18-19159
Description: K2971
Operator: Krüger / Theissmann
Submitter: Müller
Date: 25.05.2022

Preparation: polished section, uncoated
Measurement: SEM, 5kV, 7mm, ETD
Resolution: 5.8 nm / pix

Pigment Size



Median, number based:	0,230 µm	Fraction > 0.6 µm:	1%
Distribution width s:	0,17	Fraction < 0.2 µm:	36%
Mode, number based:	0,231 µm	Fraction < 0.1 µm:	1%
Median, mass based:	0,360 µm		

all number based values calculated based on the measured data; all mass based values are calculated based on fitted data; the modal value is strongly dependent on the bin-size distribution; a logarithmic binning scheme is used

Leverkusen,
25.05.2022

prepared:

approved: